

Regional Medical Supplies Catalogue List

Revised July 2024

Table of Contents

| Dressings, Bandages, Tape, Compression | 1 |
|---|----------|
| & Wraps | |
| Wound Care | 6 |
| Absorptives | 6 |
| Gelling Fiber & Antimicrobials Silver Gelling Fiber | 6 |
| Alginate & Antimicrobial Silver Alginate Dressings | 7 |
| Antimicrobial & Silver Dressings | 7 9 |
| Composite & Multifunction Dressings Clear Acrylic Dressings | 10 |
| Foams | 10 |
| Hydrocolloids | 11 |
| Hydrogel | 11 |
| Films & Membranes | 11 |
| Odour Absorbent Dressings | 12 |
| Non-adherent Offloading Devices | 12 13 |
| Other, Specialty & Kits | 14 |
| Negative Pressure Wound Therapy (NPWT) | 17 |
| | 20 |
| Cleansers & Antiseptic Solutions Infusion | _• |
| I.V. Supplies | 21 |
| Needles & Syringes | 21 |
| Prefilled Syringes and I.V. Solutions I.V. Catheters and Access | 22 |
| I.V. Gravity Sets & Extension Sets | 23 23 |
| I.V. Securement | 23 |
| Infusion Pump & Supplies | 24 |
| Elastomeric Disposable Pump & Supplies | 26 |
| Blood Collection Supplies | 26 |
| Other IV Supplies & Accessories | 27 |
| I.V. Kits | 28 |
| Ostomy Supplies and Accessories | 37 |
| Ostomy Skin Barriers for 2 Piece System | 37 |
| Urostomy Pouches for 2 Piece System | 38 |
| Drainable Pouches for 2 Piece System | 38 |
| One Piece Urostomy System | 39 |
| One Piece System with Drainable Pouches | 39 |
| Pediatric 2 Piece System | 40 |
| Ostomy Accessories | 40 |
| Urinary Supplies | 42 |
| Foley Catheters | 42 |
| Intermittent Catheters | 43 |
| Urinary Drainage Bags Catheter Stabilization | 43 44 |
| Catheter Accessories | 44 |
| Continuous Bladder Irrigation | 44 |
| Urinary Kits | 45 |
| Feeding, Suction & Tracheostomy Supplies | 47 |
| Enteral Feeding | 47 |
| Suction Catheters, Tubing and Accessories | 48 |
| Tracheostomy Products | 49 |
| Oxygen Monitoring | 50 |
| Protective Wear | 51 |

| Appendix A - Depot/Drop-off Locations | A |
|---------------------------------------|----|
| Algoma | A1 |
| Nipissing - Temiskaming | A2 |
| Cochrane | A4 |
| Sudbury-Manitoulin-Parry Sound | A6 |

- The Regional Medical Supplies Catalogue List identifies the items funded by Home and Community Care Support Services North East that are available to eligible patients who are authorized for Home and Community Care nursing services as per applicable policies.
- If different products are requested, an automatic substitution policy may be initiated to provide the equivalent product as identified within this list.
- The Home and Community Care Support Services North East is not obligated to fund or provide any items that are not identified on this list.
- Home and Community Care Support Services North East-funded nurses are expected to be trained on the supplies identified on this list, and use the supplies in alignment with the manufacturers' recommendations.

Searching the Medical Supplies Catalogue List Electronically

BOOKMARKS: Make sure the Navigation Pane is showing on the left hand side. If it is not press <F4> and then click on the second Navigation Pane option that looks like a ribbon. This will display the bookmarks. Click on the bookmark you want to view.

PRODUCTS: Click on the **Edit menu> Find>** or press **<Ctrl + F >** and type in any word pertaining to product such as product code, description, manufacturer or manufacturer code and then press Enter.



To navigate within the document, click the arrows or enter in a specific page number



Please report any errors or feedback to NE.Contracts@ontariohealthathome.ca

DRESSINGS, BANDAGES, TAPE, COMPRESSION & WRAPS

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | | |
|---------------------|--|-------------------|----------------------------|-----------------------------------|---------|------------------|--|--|--|
| materials in t | | cut and placed in | a wound because le | oose fibers (lint) r | nay ge | t lost in the w | o the fibers, wicked away, or passed through into other absorbent yound and delay healing. Nonwoven gauze consists of fibers pressed rbent capacity. | | |
| ROLL BAN | NDAGES: Knitted stretch gauze in roll configurat | ion used as seco | ndary dressing for li | mbs and to hold o | dressir | ıgs in place. I | ncludes sterile and non-sterile. | | |
| SDR-0020 | Sterile Small Woven Gauze 2" Roll | Gazetex | Derma Sciences DUP9321 | 2.25" x 3 yards 5.7 cm x 2.7 m | Υ | 7 | Gazetex bandage rolls are washed, low-linting woven gauze rolls; 6-ply. These bandages are used for wrapping or packing large and deep wounds. When used as a primary dressing, the dressing provides fast | | |
| SDR-0021 | Sterile Large Woven Gauze 4" Roll | Gazetex | Derma Sciences DUP9322 | 4.5" x 4 yards 11.4cm x 3.7m | Υ | | wicking, absorbency, and aeration. When used as a secondary dressing, the dressing provides additional cushioning and comfort. Does not contain natural rubber latex | | |
| SDR-0007 | Conforming Stretch Bandage Wrap 2" (Roll) | Duform | Derma Sciences DUP75102 | 2" x 4.1 yards 5 cm x 3.75 m | Υ | 7 | Synthetic knitted stretch bandage ideal as a secondary dressing. Non- | | |
| SDR-0008 | Conforming Stretch Bandage Wrap 4" (Roll) | Duform | Derma Sciences DUP75104 | 4"x 4.1 yards 10cm x 3.75m | Υ | 7 | Sterile | | |
| SDR-0060 | Conforming Cotton Bandage 4" (Roll) | Dutex | Derma Sciences DUP76783 | 4" x 4.1 yards 10cm x 3.75m | Υ | 7 | Cotton, 2 ply woven bandage ideal as a secondary dressing. Nonstretch and Non-sterile | | |
| DRESSING | GS AND PACKING: Sterile dressing products | to cover surgical | incisions to control I | oleeding, absorb | exudat | e, provide pr | otection, and non-sterile products to secure primary dressings in place. | | |
| SDR-0013 | Oval Eye Pad Non-woven Sterile | Surecare | Derma Sciences DUP82911 | 1-5/8" x 2-5/8" 4 x 6.5 cm | Υ | 7 | | | |
| SDR-0001 | Abdominal Pad Sterile 5"x9" | Dupad | Derma Sciences DUP87059 | 5" x 9" 12.7 x 22.8 cm | Υ | 7 | Non-woven; non-adherent | | |
| SDR-0002 | Abdominal Pad Sterile 8"x10" | Dupad | Derma Sciences DUP87810 | 8" x 10" 20 x 25 cm | Υ | 7 | Non-woven; non-adherent | | |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | Ducare | Derma Sciences DUP82212 | 2" x 2" 5 x 5 cm | Υ | 14 Packs of 2 | Cotton, 12 Ply, Pack of 2 | | |
| SDR-0004 | Gauze Sterile Woven 4"x4" (2 per Pack) | Ducare | Derma Sciences DUP82408 | 4" x 4" 10 x 10cm | Υ | 14 Packs of 2 | Cotton, 8 Ply, Pack of 2 | | |
| SDR-0017 | Sterile Packing Plain 1/4" (Bottle) | PAK-ITS | Derma Sciences DUP59120 | 1/4" x 5 yards 0.64cm x 4.5m | Υ | 2 | Single-ply ribbon gauze strips that can be used dry to absorb drainage | | |
| SDR-0018 | Sterile Packing Plain 1/2" (Bottle) | PAK-ITS | Derma Sciences DUP59220 | 1/2" x 5 yards 1.27cm x 4.5m | Υ | 2 | in exuding wounds or soaked with a medication of choice. A secondary dressing is needed to hold the strips in place and to protect the wound from outside contaminants. | | |
| SDR-0019 | Sterile Packing Plain 1" (Bottle) | PAK-ITS | Derma Sciences DUP59320 | 1" x 5 yards 2.54cm x 4.5m | Υ | 2 | from outside contaminants. (Bottles not for use at clinic - See single use packs below). | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|--|------------------|-----------------------------|---------------------------|---------|------------------|--|
| SDR-0022 | Sterile Gauze Packing Cotton 2x36" | PAK-ITS | Derma Sciences DUP564036 | 2" x 36" 5 cm x 90 cm | Υ | 3 | |
| SDR-0088 | Sterile Gauze Packing Cotton Single Use 1/4" | PAK-ITS | Derma Sciences DUP561036 | 1/4" x 36" 0.64 x 90cm | Υ | 3 | |
| SDR-0089 | Sterile Gauze Packing Cotton Single Use 1/2" | PAK-ITS | Derma Sciences DUP4853 | 1/2" x 54" 1.3 x 137cm | Υ | 3 | |
| SDR-0090 | Sterile Gauze Packing Cotton Single Use 2" | PAK-ITS | Derma Sciences DUP59420 | 2" x 54" 5 cm x 4.5m | Υ | 3 | |
| SDR-0023 | Sterile Gauze Packing Cotton Fan Folded 2x72" | PAK-ITS | Derma Sciences DUP4869 | 2" x 72" 5 cm x 182 cm | Υ | 3 | |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per Pack) | DuSoft | Derma Sciences DUP84122 | 2" x 2" 5 x 5 cm | Υ | 14 Packs of 2 | Synthetic, 4 Ply, Pack of 2 |
| SDR-0065 | Sponge Synthetic Sterile Non-Woven 4"x4" (2 per Pack) | DuSoft | Derma Sciences DUP84144 | 4" x 4" 10 x 10cm | Υ | 14 Packs of 2 | Synthetic, 4 Ply, Pack of 2 |
| SDR-0076 | Sponge Synthetic Non-sterile Non-Woven 4"x4" | DuSoft | Derma Sciences DUP94144 | 4" x 4" 10.2 x 10.2 cm | Υ | 1 Pack of 200 | Synthetic, 4 Ply, Pack of 200 |
| Available in a | | The non-woven | sponges provide su | | | | g around a tracheostomy tube or any type of tube, IV line/PICC line. fast-wicking, and non-adherent to wound beds. The woven gauze |
| SDR-0055 | Drain Gauze Sponge Tracheostomy Sterile Woven 4"x4" | Ducare | Derma Sciences DUP82910 | 4" x 4" 10 x 10cm | Υ | 14 | Woven, 12 Ply, Pack of 1 |
| SDR-0058 | Drain Sponge Sterile Non-woven 4x4 | DuSoft | Derma Sciences DUP84916 | 4" x 4" 10 x 10cm | Υ | 14 Packs of 2 | 6 Ply Non-woven, Pack of 2 |
| SOT-0036 | Drain Sponge Sterile Non-woven 2x2 | DuSoft | Derma Sciences DUP84918 | 2" x 2" 5 x 5 cm | Υ | 14 Packs of 2 | 6 Ply Non-woven, Pack of 2 |
| SOT-0024 | Drain Sponge Antimicrobial Sterile Non-woven 4x4 | Excilon | Covidien Z7088 | 4" x 4" 10 x 10cm | Υ | 14 Packs of 2 | PHMB, 6 Ply Non-woven, Pack of 2 |
| TAPES & S | SECUREMENT: Medical and surgical tapes use | ed as a secureme | ent aid in dressing o | r tubings (NG, IV | , Foley | , etc.). | |
| SDR-0039 | Non Allergenic Self-Adhesive Dressing Retention Tape 2" | Ultrafix | Derma Sciences DUP62036 | 2" x 1 yrd 5 cm x 1 m | Υ | 7 m | A non-sterile and non-woven fixative tape that is breathable and gentle to the skin. Custom cut by metre lengths |
| SDR-0040 | Non Allergenic Self-Adhesive Dressing Retention Tape 4" | Ultrafix | Derma Sciences DUP62037 | 4" x 1 yrd 10 cm x 1 m | Υ | 7 m | A non-sterile and non-woven fixative tape that is breathable and gentle to the skin. Custom cut by metre lengths |
| SDR-0048 | Montgomery Straps and Ties | Montgomery | Derma Sciences DUP97991 | 18 cm x 28 cm | Υ | 1 | Non-sterile, non-woven |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|---|---------------------------|-------------------------------|---------------------------------|-----|---------------|---|
| SDR-0086 | White Surgical Tape 1" Single Use Roll | Transpore (or equivalent) | 3M 1534S-1 (or equivalent) | 1" x 1.5 yard 2.5 cm x 1.3 m | Υ | 1 | Universal hypoallergenic tape can be used for all applications. |
| SDR-0087 | Blue Silicone Tape 1" Single Use Roll | Micropore S | 3M 2770S-1 (or equivalent) | 1" x 54" 2.5 cm x 1.3 m | Υ | | Only order for patients with very fragile at risk skin who react to SDR-0086 1" White Surgical Tape |
| SDR-0037 | Wound Closure Strip 1/2" (6 Strips) Sterile | Primaflex | priMED Medical PM3-71241 | 1/2" x 4" 12 x 100 mm | Υ | 2 Packs | 6 Strips per pack |
| SDR-0038 | Wound Closure Strip 1/4" (3 Strips) Sterile | Primaflex | priMED Medical PM3-71431 | 1/4" x 3" 6 x 75 mm | Υ | 2 Packs | 3 Strips per pack |

BANDAGES & WRAPS: Materials that are dry woven or non-woven and may stretch and conform to body contours. Fabric composition may include cotton, polyester, rayon or nylon. They function for cleansing, packing or to provide absorption as a second layer or dressing, to hold a wound cover in place, to apply pressure or to cushion a wound. Some are sterile or non-sterile in bulk and may be adhesive. Other uses include securing IVs and splints and providing mild compression/support for strains, sprains and edema.

| SDR-0027 | Tubular Net Elastic Bandage - Fingers/Toes/Small Wrist | Medline Elastic Net | NONNET02 | 7 2/3" 19.5 cm | Υ | 2 x 1 metre | |
|----------|---|------------------------|----------|--------------------|---|-------------|--|
| SDR-0083 | Tubular Net Elastic Bandage - Wrist & Elbow | Medline Elastic Net | NONNET03 | 9 1/3" 23.7 cm | Υ | 2 x 1 metre | |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | Medline Elastic Net | NONNET04 | 11 2/3" 29.6 cm | Υ | 2 x 1 metre | |
| SDR-0029 | Tubular Net Elastic Bandage - Knee & Thigh | Medline Elastic Net | NONNET05 | 12 3/8" 31.4 cm | Υ | 2 x 1 metre | Non-adhesive and holds dressing/devices securely in place while allowing freedom of movement. Roll of 25 yards cut to 1 meter lengths |
| SDR-0084 | Tubular Net Elastic Bandage - Shoulder Hip Thigh | Medline Elastic Net | NONNET06 | 17" 61 cm | Υ | 2 x 1 metre | , , |
| SDR-0030 | Tubular Net Elastic Bandage - Chest Back Perineum Axilla | Medline Elastic Net | NONNET07 | 28 2/3" 72.8 cm | Υ | 2 x 1 metre | |
| SDR-0085 | Tubular Net Elastic Bandage - Chest Region Large | Medline Elastic Net | NONNET08 | 31 1/3" 79.6 cm | Υ | 2 x 1 metre | |
| SDR-0031 | Tubular Net Elastic Bandage - Chest Region XL | Medline Elastic Net | NONNET09 | 35" 88.9 cm | Υ | 2 x 1 metre | Non-adhesive and holds dressing/devices securely in place while allowing freedom of movement. SPECIAL ORDER ITEM - May take up |
| SDR-0035 | Tubular Net Elastic Bandage - Chest Region XXL | Medline Elastic Net | NONNET10 | 41" 93.1 cm | Υ | 2 x 1 metre | to 5 days for stock. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Max | Notes |
|---------------------|--|-------------------|------------------------------|----------------------|--------|---------------|---|
| SCM-0003 | Mesh Pants XL / XXL | Wings | Cardinal Health Z707A | 40"-56" | Υ | 1 | Incontinuos augustica are not funded as nor policy. These mach ponto |
| SCM-0004 | Mesh Pants Large / XL | Wings | Cardinal Health Z706A | 32"-42" | Υ | 1 | Incontinence supplies are not funded as per policy. These mesh pants are only authorized for patients requiring this product as a means of containing wound care dressings to coccyx, perineal area, etc. |
| SCM-0005 | Mesh Pants Small / Medium | Wings | Cardinal Health Z705A | 22"-32" | Υ | 1 | - bornairing wearing dare drobbings to boodyx, perinted area, etc. |
| SDR-0068 | Multipurpose Tubular Elastic Bandage Size C | Reliamed | Cardinal Health ZG275TB | 2.75" x 36" | Υ | 2 | - Limb sizes C through G available. Other sizes may available only as an exception item special order with approval. Sizes intended as a |
| SDR-0069 | Multipurpose Tubular Elastic Bandage Size D | Reliamed | Cardinal Health ZG30TB | 3" x 36" | Υ | 2 | guideline only and should be used in conjunction with existing clinical protocols that consider physiological state and wound and skin |
| SDR-0070 | Multipurpose Tubular Elastic Bandage Size E | Reliamed | Cardinal Health ZG35TB | 3.5" x 36" | Υ | 2 | condition Order in 1 yard lengths (36") |
| SDR-0071 | Multipurpose Tubular Elastic Bandage Size F | Reliamed | Cardinal Health ZG40TB | 4" x 36" | Υ | 2 | - If overlapped, applies pressure of 10-15 mm Hg - Can be washed (gentle cycle & air dry) and reused multiple times |
| SDR-0072 | Multipurpose Tubular Elastic Bandage Size G | Reliamed | Cardinal Health ZG45TB | 4.5" x 36" | Υ | 2 | - Contains latex |
| SVS-0003 | High Stretch Elastic Bandage Non-Latex 4" | Dusor | Derma Sciences DUP71414LF | 4" x 5 yards | Υ | 1* | *May order 2 for initial order only if needed. Non-Sterile |
| SVS-0004 | High Stretch Elastic Bandage Non Latex 6" | Dusor | Derma Sciences DUP71416LF | 6" x 5 yards | Υ | 1* | *May order 2 for initial order only if needed. Non-Sterile |
| COMPRES | SION WRAPS: Short Stretch multi-layer elastic | c stretch bandage | e system with limited | d extensibility (les | s than | 100%) for u | se to provide moderate to high compression |
| SVS-0006 | 2 Layer Compression System 35-40mmHg | Coban 2 | 3M 2094 | Each | Υ | 2 | 2 Rolls - Comfort layer and compression layer. Ideal for majority of patients with venous leg ulcers, lymphedema and other conditions where compression bandage therapy is appropriate. |
| SVS-0007 | 2 Layer Compression System 25-30 mmHg | Coban 2 Lite | 3M 2794 | Each | Υ | 2 | 2 Rolls - Comfort layer and compression layer. Designed to be more comfortable for patients less tolerant of compression bandages. |
| SVS-0009 | Tubular Elastic Compression Small Navy Stripe | EdemaWear | CA600001 | Each Pair (2) | Υ | 1/Month | Up to 2 Pair at admission then 1 Pair per month as needed |
| SVS-0010 | Tubular Elastic Compression Medium Yellow Stripe | EdemaWear | CA960001 | Each (1) | Y | 1/Month | 15-20 mmHg; Air-permeable; Can be cut Small (Navy): Up to 18"/45cm circumference (Foot/Calf/Arm) Medium (Yellow): Up to 30"/75cm circumference (Knee/Thigh/Arm) |
| SVS-0011 | Tubular Elastic Compression Large Red Stripe | EdemaWear | B120L01 | Each (1) | Υ | 1/Month | Large (Red): Up to 46"/115cm circumference (Thigh/Torso) |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | | |
|---------------------|---|-----------------|-----------------------|-----------|-----|---------------|---|--|--|
| OTHER & | SPECIALTY: | | | | | | | | |
| SCM-0019 | Blue Absorbent Pad 23"x24" | Cardinal Health | Cardinal UPSMD2324 | 23" x 24" | Y | 1/Visit | Incontinance supplies are not funded as per approved policy. This product is authorized only for patients for containment of drainage or for methotrexate administration in-home/clinic ONLY. Approved clinics may order in bulk. | | |
| DRESSING | G KITS: | | | | | | | | |
| SDR-0005 | Dressing Tray Kit - Disposable Single Use | MED-Rx | CHS 85-4009 | Each Kit | Y | 1 | For catheter treatments only (PNT or PICC). Not for use in wound care. | | |
| 3 Compa | 3 Compartment 1200 CC Outer Tray | | | | | | | | |
| Rayon C | otton Balls Large Absorbent | | | 6 | | | | | |
| Gauze 4 | " x 4" - 8 Ply | | | 4 | | | | | |
| 4" Plastic | c Transfer Forcep | | | 3 | | | | | |
| Absorbe | Absorbent Towel - 3 Ply | | | | | | | | |
| CSR Ste | CSR Sterile Field Towel | | | | | | | | |
| Waste B | ag with Tie | | | 1 | | | | | |

WOUND CARE

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|---|----------------|----------------------|-----------------|-------|---------------|--|
| Oraci ocac | | | Oddc | | | IVIUA | |
| Home and C | Communuty Care Support Services North E | ast-funded pro | ofessional nurse | es are required | to be | trained of | on all wound care supplies that are funded. Nurses are expected to |

Home and Communuty Care Support Services North East-funded professional nurses are required to be trained on all wound care supplies that are funded. Nurses are expected to use the products in alignment with defined Wound Pathways and the manufacturer's product recommendations for use, contraindications, wear time, etc.

ABSORPTIVES: Multi-layer wound covers that provide either a semi-adherent quality or a non-adherent layer, combined with highly absorptive layers of fibers, such as cellulose, cotton or rayon. Designed to minimize adherence to the wound and manage exudate. They may be used as a primary or secondary dressing to manage surgical incisions, lacerations, abrasions, burns, donor or skin graft sites, or any exudating wound.

| SWC-0107 | Absorbent Dressing 4"x6" | Mesorb | Mölnlycke 677201 | 4" x 6" 10 x 15 cm | Υ | 7 | | |
|----------|---------------------------------|----------------------------|---------------------|------------------------|---|---|--|--|
| SWC-0108 | Absorbent Dressing 6"x8" | Mesorb | Mölnlycke 677601 | 6" x 8" 15 x 20 cm | Υ | 7 | | |
| SWC-0109 | Absorbent Dressing 8"x12" | Mesorb | Mölnlycke 677801 | 8" x 12" 20 x 30 cm | Υ | 7 | | |
| SWC-0142 | Super-Absorbent Dressing 4"x6" | Mextra Super- absorbent | Mölnlycke 610100 | 4" x 6" 10 x 15 cm | Υ | 4 | - For best result, the pad should overlap the dry | |
| SWC-0143 | Super-Absorbent Dressing 6"x8" | Mextra Super- absorbent | Mölnlycke 610300 | 6" x 8" 15 x 20 cm | Υ | 4 | surrounding skin by at least 2 cm. - May be left in place for several days. The dressing should be changed according to the clinical condition of the wound, or when saturated. | |
| SWC-0144 | Super-Absorbent Dressing 8"x12" | Mextra Super- absorbent | Mölnlycke 610500 | 8" x 12" 20 x 30 cm | Υ | 4 | | |

GELLING FIBER & ANTIMICROBIAL SILVER GELLING FIBER DRESSINGS: The products have a "primary mode of action" to absorb wound exudate and maintain a moist wound environment • Converts to a solid gel when activated by moisture • Includes synthetic gelling fiber dressings with and without silver antimicrobial properties • Supports autolytic debridement • Best for moderate amounts of exudate • Should not be used on dry wounds.

| SWC-0150 | Gelling Fiber Dressing 4"x4" | Aquacel Extra | ConvaTec 420672 | 4" x 4" 10 x 10 cm | Υ | 4 | | |
|----------|--|----------------------|--------------------|------------------------|---|---|--------------|----|
| SWC-0151 | Gelling Fiber Dressing 6"x6" | Aquacel Extra | ConvaTec 420673 | 6" x 6" 15 x 15 cm | Υ | 4 | | |
| SWC-0106 | Gelling Fiber Ribbon 1"x18" | Aquacel | ConvaTec 403770 | 1" x 18" 2 x 45 cm | Υ | 4 | | |
| SWC-0152 | Gelling Fiber Antimicrobial Silver Dressing 4"x4" | Aquacel AG+ Extra | ConvaTec 413567 | 4" x 4" 10 x 10 cm | Υ | 4 | | FI |
| SWC-0153 | Gelling Fiber Antimicrobial Silver Dressing 8"x12" | Aquacel AG+ Extra | ConvaTec 413569 | 8" x 12" 19 x 29 cm | Υ | 4 | | |
| SWC-0032 | Gelling Fiber Antimicrobial Silver Ribbon 1"x18" | Aquacel AG+ | ConvaTec 413571 | 1" x 18" 2 x 45 cm | Υ | 4 | Ionic Silver | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|--------------------|----------------------------|------------------------------|-----|---------------|---|--|
| exudate, such | | ers and venous in | sufficiency ulcers • | Has hemostatic | | | nposed of seaweed derivatives • Indicated for wounds wures a secondary dressing • Should not be used on dr | |
| SWC-0036 | Alginate Dressing 4"x4" | Tegaderm | 3M 90112 | 4" x 4" 10 x 10 cm | Υ | 3 | | |
| SWC-0120 | Antimicrobial Silver Alginate Dressing | Silvercel | KCI Medical CAD011 | 4.3" x 4.3" 11 x 11 cm | Υ | 3 | | |
| SWC-0154 | Antimicrobial Silver Rope Alginate Dressing | Silvercel | KCI Medical CAD7230 | 2.5 x 30.5cm | Υ | 3 | | |
| and surgical i | | | | | | | risk of infection in partial- and full-thickness wounds, or product. Does not replace systemic antibiotics for deep | |
| SWC-0034 | Non-adherent Antiseptic Impregnated Gauze 4"x4" | Bactigras | Smith & Nephew 66003650 | 4" x 4" 10 x 10 cm | Y | 4 | For adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers | |
| SSP-0003 | Ointment with lodine 10Gr Tube | lodasorb | Smith & Nephew 66060630 | 10 Gram Tube | Υ | 1 | Cadexomer lodine should not be used in patients with known or suspected lodine sensitivity, history of any thyroid disorder, or pregnant or lactating women | IODOSOR® LOCACCETE ISFA Commit Copyritor Iodosomber grade Copyritor to downship grade |
| SDR-0078 | Non-adhesive Povidone-lodine Dressing 2"x2" | Inadine | KCI Medical P01481 | 2" x 2" 5 x 5 cm | Υ | 3 | 10% Povidone lodine; equivalent to 1.0% available iodine | |
| SDR-0079 | Non-adhesive Povidone-Iodine Dressing 3.75"x3.75" | Inadine | KCI Medical P01512 | 3.75 x 3.75" 9.5 x 9.5 cm | Υ | 3 | 10% Povidone lodine; equivalent to 1.0% available iodine | |
| SWC-0119 | Silver Antimicrobial Barrier Dressing Non- Absorbent 4"x4" | Acticoat Flex 3 | Smith & Nephew 66800399 | 4" x 4" 10 x 10cm | Υ | 3 | Not compatible with oil-based products, such as petrolatum | |
| SWC-0156 | Antimicrobial Silver Mesh | Acticoat Flex 7 | Smith & Nephew 66800397 | 4" x 5" 10 x 12.5cm | Υ | 1 | Up to 7 day wear time. Not to be ordered with dressing changes more frequently than Q3 days. Can be cut to fit. | |

dressing changes more frequently than Q3 days. Can be cut to fit.

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|-------------------------|--|----------------------------------|-----|---------------|--|--|
| SWC-0010 | PHMB Antimicrobial Gauze Square 4"x4" | Curity | Cardinal Health Z2533 | 4" x 4" 10.2 x 10.2 cm | Υ | 7 | Impregnated with 0.2% PHMB; 12 Ply | |
| SWC-0112 | PHMB Antimicrobial Non-Woven Sponges 4"x4" | Kendall | Cardinal Health Z8044AMD or equivalent | 4" x 4" 10.1 x 10.1 cm | Y | 7 | Impregnated with 0.2% PHMB | |
| SWC-0113 | PHMB Antimicrobial Gauze Roll 4.5" | Kerlix | Cardinal Health Z3332 | 4.5" x 4 1/8 yd 11.4cm x 3.7m | Υ | 7 | Impregnated with 0.2% PHMB; 6 Ply | |
| SWC-0115 | Antimicrobial Gauze PHMB Ribbon 1/4" x 1 yard | Curity | Cardinal Health Z7831AMD | 1/4" x 36" 0.6 x 91.4 cm | Υ | 7 | Impregnated with 0.2% PHMB | |
| SWC-0116 | Antimicrobial Gauze PHMB Ribbon 1/2" x 1 yard | Curity | Cardinal Health Z7832AMD | 1/2" x 36" 1.3 x 91.4 cm | Υ | 7 | Impregnated with 0.2% PHMB | 11/18 |
| SWC-0117 | Antimicrobial Gauze PHMB Ribbon 1" x 1 yard | Curity | Cardinal Health Z7833AMD | 1" x 36" 2.5 x 91.4 cm | Υ | 7 | Impregnated with 0.2% PHMB | |
| SWC-0118 | Antimicrobial Gentian Violet Methylene Blue Dressing 4x4" | Hydrofera Blue | Hydrofera HB4414 | 4" x 4" 10 x 10 cm | Υ | 3 | | |
| SWC-0137 | Antimicrobial Gentian Violet Methylene Blue Dressing 6x6" | Hydrofera Blue | Hydrofera HB6614 | 6" x 6" 15 x 15 cm | Υ | 3 | - Not to be used longer than 6 months without consuting Wound Specialty Nurse. | |
| SWC-0155 | Antimicrobial Gentian Violet Methylene Blue Moist Foam Dressing 4x5" | Hydrofera Blue Ready | Hydrofera HBRS4520 | 4" x 5" 10 x 12.5cm | Y | 1 | - Not to be used on pregnant or nursing patients. | |
| SWC-0138 | Antimicrobial Gentian Violet Methylene Blue Tunneling Dressing 6x6" | Hydrofera Blue | Hydrofera HBT0906 | 9 mm (Pack of 2) | Υ | 3 | | |
| SWC-0131 | Antimicrobial Gel with Silver 0.25 oz | SilvaSorb | Medline MSC93025EP | 0.25 oz Tube | Υ | As required | For clinic patients only | SILVASORB' GEL |
| SWC-0132 | Antimicrobial Gel with Silver 1.5 oz | SilvaSorb | Medline MSC9301EP | 1.5 oz Tube | Υ | 1 | Community patients - not clinic | - 100 2 100 41 1 |
| SWC-0134 | Textile Antimicrobial Silver (Custom order by 12" lengths) | InterDry AG | Coloplast 7910 | 10" x Custom 12" lengths | Y | 4 | Moisture-wicking fabric. Custom order per 12" lengths | InterDry More than the control reason of th |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | | |
|---------------------|--|-----------------------|----------------------|---------------------------------|-----|---------------|---|--|--|
| SWC-0157 | Foam Silver Dressing | Biatain AG | Coloplast 9622 | 4" x 4" 10 x 10 cm | Y | 4 | Non-adhesive | Research And Bustains And Busta | |
| SWC-0145 | Silver Non-Adhesive Pad Dressing | PolyMem Silver | 1044 | 4.25" X 4.25" 10.8 x 10.8 cm | Υ | 2 | Cleanses, fills, absorbs and moistens wounds. Combined primary/secondary dressing or as a secondary dressing. Can be cut to size. | | |
| SWC-0146 | Silver Non-Adhesive Pad Dressing (MAX) | PolyMem MAX Silver | 1045 | 4" X 4" 10 x 10 cm | Υ | | MAX is a thicker PolyMem configuration for wounds with heavier drainage or when longer wear time is desired. | | |

COMPOSITE & MULTIFUNCION DRESSINGS: Multilayered, combination dressings to provide multiple functions, such as a bacterial barrier, absorption and adhesion. Usually, they are comprised of multiple layers and incorporate a semi- or non-adherent pad that covers the wound. Some are appropriate for autolysis. May also include an adhesive border of non-woven fabric tape or transparent film. Composite dressings can function as either a primary or a secondary dressing on a wide variety of wounds.

| SWC-0147 | Non-Adhesive Pad Dressing | PolyMem | 5044 | 4" x 4" 10 x 10 cm | Y | 2 | Cleanses, fills, absorbs and moistens wounds. Combined primary/secondary dressing or as a secondary dressing. Can be cut to size. | |
|----------|---|---------------------|---------------------------|---------------------------------|---|---|--|--|
| SWC-0148 | StatLock Universal Plus Catheter Stabilization Device Small 6 - 8.5 FR | PolyMem MAX | 5045 | 4.25" x 4.25" 10.8 x 10.8 cm | Υ | 2 | - MAX is a thicker PolyMem configuration for wounds with heavier drainage or when longer wear time is desired. | |
| SWC-0090 | Protease Inhibitor Dressing with Collagen | Promogran | KCI M772028 | 28 cm² Hexagonal | Y | 2 | Review manufacturer's information regarding contraindications. | |
| SWC-0158 | Protease Inhibitor Dressing with Silver | Promogran Prisma | KCI PS2028 | 28 cm² Hexagonal | Y | 2 | Protease inhibitor should be discontinued if there is no progression towards wound healing (less than 10% healing in 3 weeks). | |
| SVS-0008 | Non-Elastic Zinc Paste Bandage | Primer | Derma Sciences GL300-1 | 7.6 cm x 9 m | Y | 2 | A gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide. Ideal for managing: venous ulcers, lymphaticedema, minor sprains, strains, fractures and thrombophlebitis. Contains no gelatin which prevents the hardening or caking of the bandage. Made without preservatives and with non-raveling thread locked gauze. | |

| Γ | CHRIS | Doscription | Brand Name | Manufacturer | Size | APO | Weekly | Notes |
|---|------------|-------------|------------|--------------|------|-----|--------|-------|
| | Order Code | Description | Dianu Name | Code | Size | AFU | Max | Notes |

CLEAR ACRYLIC DRESSINGS: Transparent film contact layer and clear, acrylic polymer pad, topped with breathable, waterproof film. Enables clinicians to monitor small to moderately exuding wounds without changing the dressing • Supports autolytic debridement • Extended wear time • Skin graft donor sites • Skin tears and pre-tibial lacerations • Post-operative wounds • First (superficial thickness) and second (partial thickness) degree burns and abrasions • Do not cut acrylic pad

| SWC-0139 | Absorbent Clear Acrylic Small Oval 7.6x9.5cm | Tegaderm | 3M 90800 | 3 in. x 3-3/4 in. 7.6 x 9.5 cm | Υ | For pressure ulcers, skin tears and other low to moderately exuding wounds. Dressing should | |
|----------|--|----------|-------------|-------------------------------------|---|---|--|
| SWC-0140 | Absorbent Clear Acrylic Medium Oval 11x12cm | Tegaderm | 3M 90801 | 4-3/8 in. x 5 in. 11.1 x 12.7 cm | Y | overlap the wound margins by at least 3/8" (1 cm). | |

FOAMS: Foam dressings are sheets and other shapes of foamed polymer solutions (most commonly polyurethane) with small, open cells capable of holding fluids. They may be impregnated or layered in combination with other materials. Absorption capability depends on thickness and composition. The area in contact with the wound surface is non-adhesive for easy removal. Available with an adhesive border and/or a transparent film coating that acts as a bacterial barrier. Indicated for partial- and full-thickness wounds. Used on lightly- to heavily-exudating wounds, depending on product. Occlusive foams without silver should not be used on infected wounds.

| SWC-0099 | Foam Non-Border 4"x4" | Mepilex XT | Mölnlycke 211100 | 4" x 4" 10 x 10 cm | Υ | 4 | |
|----------|--|-----------------------------|---------------------|-------------------------|---|---|---------|
| SWC-0100 | Foam Non-Border 8"x8" | Mepilex XT | Mölnlycke 211400 | 8" x 8" 20 x 20 cm | Υ | 4 | |
| SWC-0101 | Foam Non-Border Soft Silicone 6"x6" | Mepilex Lite | Mölnlycke 284300 | 6" x 6" 15 x 15 cm | Υ | 4 | |
| SWC-0102 | Foam Border Silicone 3"x3" | Mepilex Border Flex | Mölnlycke 595211 | 3 "x 3" 7.5 x 7.5 cm | Υ | 4 | |
| SWC-0103 | Foam Border Silicone 4"x4" | Mepilex Border- Flex | Mölnlycke 595311 | 4" x 4" 10 x 10 cm | Υ | 4 | |
| SWC-0104 | Foam Border Silicone 6"x6" | Mepilex Border- Flex | Mölnlycke 595411 | 6" x 6" 15 x 15 cm | Υ | 4 | |
| SWC-0105 | Foam Border Silicone 4"x4" | Mepilex Border Flex Lite | Mölnlycke 581300 | 4" x 4" 10 x 10 cm | Υ | 4 | SafetaC |
| SWC-0129 | Foam Non-Adhesive with Ibuprofen 4"x4" | Biatain IBU | Coloplast 4110 | 4" x 4" 10 x 10 cm | Υ | 4 | |
| SWC-0130 | Foam Non-Adhesive with Ibuprofen 6"x6" | Biatain IBU | Coloplast 4115 | 6" x 6" 15 x 15 cm | Υ | 4 | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|----------------------------------|---|---------------------------------------|--|--|---------------------|--------------------------|--|
| HYDROCO invasion. Abs | DLLOIDS: Wafers, powders or pastes composed sorption capability depends on thickness and comp | osition. Useful o | or carboxymethylcen areas that require | contouring, such | as hee | ve dressir els and sa | Ing, contributes to autolytic debridement. Creates occlusive barrier against bacteria cral ulcers. Indicated for partial and full thickness wounds with or without necrotice skin and should not be used on heavily draining or infected wounds. |
| SWC-0051 | Hydrocolloid Regular Square 4x4" | Comfeel Plus | Coloplast 33110 | 4" x 4" 10 x 10 cm | Υ | 3 | |
| SWC-0052 | Hydrocolloid Regular Square 6x6" | Comfeel Plus | Coloplast 33115 | 6" x 6" 15 x 15 cm | Υ | 3 | |
| SWC-0080 | Hydrocolloid Regular 8"x8" | Comfeel Plus | Coloplast 33120 | 8" x 8" 20 x 20 cm | Υ | 3 | 6 mmg |
| SWC-0056 | Hydrocolloid Regular Sacral | Comfeel Plus | Coloplast 33285 | 7" x 8" 18 x 20 cm | Υ | 3 | |
| SWC-0045 | Hydrocolloid Transparent Square 4x4" | Comfeel Plus Transparent | Coloplast 33533 | 4" x 4" 10 x 10 cm | Υ | 3 | |
| SWC-0048 | Hydrocolloid Transparent Square 6x6" | Comfeel Plus Transparent | Coloplast 33539 | 6" x 6" 15 x 15 cm | Υ | 3 | |
| SWC-0047 | Hydrocolloid Transparent Rectangle 3.5x5.5" | Comfeel Plus Transparent | Coloplast 33536 | 3.5" x 5.5" 9 cm x 14 cm | Υ | 3 | |
| SWC-0053 | Hydrocolloid Transparent Oval | Comfeel Plus Transparent | Coloplast 33283 | 3.5" x 4" 9 x 11 cm | Υ | 3 | |
| rehydrate woi | | dded into gauze. | Can decrease pair | | | | d maintain a moist healing environment. The high moisture content serves to kness wounds, wounds with necrosis, minor burns and radiation tissue damage. |
| SWC-0023 | Hydrogel (Low Pectin, Low Saline) | Intrasite Gel | Smith & Nephew 7311 | 15 g | Y | 3 | |
| moisture vapo Indicated for p | our and atmospheric gases. Moisture vapour trans | mission rate (MV wounds with necre | TR) varies from film osis, and as both a | n to film. Can help primary or second | o reduc dary dre | e friction t | Ihesive. They are impermeable to liquid, water and bacteria but permeable to to susceptible skin (e.g. heels). Transparency allows visualization of the wound. so used to cover IV sites, donor sites, lacerations, abrasions and second-degree |
| SWC-0072 | Transparent Film Dressing 10x12cm | Tegaderm | 3M 1626W | 4" x 4.75" 10 x 12 cm | Υ | 3 | |
| SWC-0128 | Transparent Film Dressing 15x20cm | Tegaderm | 3M 1628 | 6" x 8" 15 x 20 cm | Υ | 3 | |
| SWC-0159 | Transparent Film Dressing 10cm x 1m | Tegaderm | 3M 16004 | 10 cm x 1 m | Υ | 1 Meter | Roll cut in 1 meter lengths per patient order |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Max | Notes | |
|---------------------|--|-----------------------------|--------------------------|------------------------------|---------|---------------|--|----------------------|
| ODOUR A | BSORBENT DRESSINGS: Odour absorbent | charcoal contain | ed within product th | nat masks the odd | our. En | sure that | dressing edges are sealed to control odour and treat th | e cause if possible. |
| SWC-0037 | Activated Charcoal Dressing with Silver 4"x4" | Actisorb Plus Silver 220 | KCI Systagenix MAS105 | 4" x 4" 10.5 x 10.5 cm | Υ | 2 | Silver is not released into the wound, therefore can | 43 |
| SWC-0135 | Activated Charcoal Dressing with Silver 2.5 x 3.75" | Actisorb Plus Silver 220 | KCI Systagenix MAS065 | 2.5" x 3.75" 6.5 x 9.5 cm | Υ | 2 | be use on silver-sensitive patients. Remains effective to reduce odour when wet or dry and can be left in place for up to 7 days. NOTE: This product | |
| SWC-0136 | Activated Charcoal Dressing with Silver 4" x 7.5" | Actisorb Plus Silver 220 | KCI Systagenix MAS190 | 4" x 7.5" 10.5 x 19 cm | Υ | 2 | should not be cut | |
| | ERENT: Porous sheets of dressings with low adh plication of topical preparations. Use with wounds | | | | | | exudate to secondary dressing. May be composed of | silicone or tulles. |
| SWC-0003 | Non Adherent Sterile Impregnated Gauze 3"x3" | Adaptic | Systagenix 2012 | 3" x 3" 7.6 x 7.6 cm | Y | 4 | | |
| SWC-0076 | Non Adherent Sterile Impregnated Gauze 3"x8" | Adaptic | Systagenix 2015 | 3" x 8" 7.6 x 20.3 cm | Υ | 7 | | |
| SWC-0160 | Non Adherent Sterile Silicone Dressing 3"x4" | Mepitel | Mölnlycke 290710 | 3"x 4" 7.5 x 10 cm | Υ | 3 | | |
| SWC-0086 | Non Adherent Sterile Impregnated Gauze 4"x4" | Jelonet | Smith & Nephew 7404 | 4"x4" 10 x 10 cm | Υ | 4 | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Max | Notes |
|---------------------|--|------------|------------------------|--|--------|---------------|--|
| | ING DEVICES: Offloading devices are used in a superior authorized prescriber (Nurse Practicioner, NSWC | | | | Contac | t Casts (| TCC), custom braces and ankle and/or foot orthoses. These products can b |
| SWC-0165 | Diabetic Walker Boot System - Small | AirCast XP | DJO Canada 01PD-S | Shoe Size (S) Male 4-7 Female 5-8 | Y | 1 | - The XP Diabetic Walker Boot System incorporates |
| SWC-0164 | Diabetic Walker Boot System - Medium | AirCast XP | DJO Canada 01PD-M | Shoe Size (M): Male 7-10 Female 8-11 | Y | 1 | an optimized rocker sole, semi-rigid shells, an innovative dual density insole and four aircells for maximum off -loading and "total contact" fit of the diabetic foot. Two extra insoles (including the dual- |
| SWC-0163 | Diabetic Walker Boot System - Large | AirCast XP | DJO Canada 01PD-L | Shoe Size (L): Male 10-13 Female 11-15 | Y | 1 | laminate Impax™ Grid insoles (including the dual- laminate Impax™ Grid insole), three stockings and one hand bulb with pressure gauge are included. - XL available only by special order and may take 3-5 |
| SWC-0016 | Diabetic Walker Boot System - X Large | AirCast XP | DJO Canada 01PD-XL | Shoe Size (XL) Male 13+ Female 15+ | Y | 1 | days for stock |
| SWC-0017 | Diabetic Walker Insole Replacement Small | AirCast XP | DJO Canada 01PI-S | Small | Υ | 1 | The initial order of the Diabetic Walker Boot System |
| SWC-0027 | Diabetic Walker Insole Replacement Medium | AirCast XP | DJO Canada 01PI-M | Medium | Υ | 1 | includes two extra insoles. For use with the AirCast XP Diabetic Walker Boot |
| SWC-0041 | Diabetic Walker Insole Replacement Large | AirCast XP | DJO Canada 01PI-L | Large | Υ | 1 | XL size available only by special order and may take |
| SWC-0049 | Diabetic Walker Insole Replacement X-Large | AirCast XP | DJO Canada 01PI-XL | Extra Large | Υ | 1 | 3-5 days for stock |
| SWC-0162 | Replacement Inflation Device for Diabetic Walker Boot | AirCast XP | DJO Canada 01PG | Each | Y | 1 | Replacement hand bulb with pressure gauge. A single hand bulb with pressure gauge is already included with each order of a Diabetic Walker Boot System. |
| SWC-0166 | Total Contact Cast Kit | Cutimed | BSN Medical 7349900 | 1 Kit | Υ | 1 | Kit includes: Stockinette x1, Adehesive Felt x1, Adhesive Foam x1, Padding Roll x2, Soft Fiberglass Cast Tape Rolls x2, Plaster of Paris Slab x1, Fiberglass Cast Tape Roll x3 |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|------------------------|------------------------|-----------------------|-----|----------------------------|--|--|
| OTHER & S | SPECIALTY: | | | | | | | |
| SWC-0161 | Hydrophilic Wound Paste | Triad | Coloplast 1964 | 71 g Tube | Υ | 1 Tube every 2 weeks | Primary mode of action to maintain a moist wound environment, wound debridement and absorbing wound exudate. Sterile and composed zinc oxide paste | The production of the producti |
| SDR-0011 | Sterile Single Use Scissors | MED-Rx | CHS 85-4202 | 4 1/2" | Υ | 7 | Sharp/Blunt. For single-use, sterile techniques only. Not to be reused. Clinic may order 1 / Patient / Visit In-Home may order 1 / Treatment Visit | 20 |
| SDR-0036 | Non-Sterile Reusable Scissors | Fisher Scientific | ALMP12 | 4 1/2" | Y | 1 per Patient Admit | Non-sterile, Single patient use. For non-invasive, non-critial use. May be re-used by the same patient after cleaning or low-level disinfection as per nursing guidelines. Not for use in nursing clinics - Use SDR-0011 Sterile Single Use Scissors instead | |
| SDR-0012 | Sterile Plastic Thumb Forceps | MED-Rx | CHS 85-4200 | Each | Υ | 1 / Visit | Sterile Single Use only 1 / Patient / Visit | |
| SDR-0025 | Skin Staple Remover | Medline Canada | Medline DYNJ04058 | Each | Υ | 1 | Single use; Sterile; Tweezer Style | |
| SDR-0062 | Wound irrigation tip for syringe | MED-RX | MED-RX 20-1104 | 1 Tip | Υ | 7 | Single use. Syringe must be ordered separately | |
| SDR-0080 | Tissue Forceps Metal Sterile Disposable | MED-Rx | CHS 85-9000 | 1 x 2 teeth 4-3/4" | Y | 1/Visit | For sharp/surgical debridement and should only be performed by clinicians with this expertise. | |
| SDR-0081 | Curved Iris Scissors Sterile | MED-Rx | CHS 85-9411 | 4 1/2" | Υ | 1/Visit | For sharp/surgical debridement and should only be performed by clinicians with this expertise | 90 |
| SWC-0059 | Wound Culture Swab | BD BBL Culture Swab | BD 220105 | Each | Υ | 1 | | COS ES Common CCC 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 |
| SWC-0064 | Single Use Safety Scalpel - Size 10 | Miltex | Integra 4-510 | Size 10 | Υ | 3 | Sharp/surgical debridement should only be performed by clinicians with this expertise. | |
| SWC-0066 | Flexible Silver Nitrate Caustic 6" Applicator Stick (3 Pack) | AMG | AMG Medical 118-395 | Each | Υ | 1 pkg 3 | Silver Nitrate 75%, Potassium Nitrate 25% | Silver Notrale Application———————————————————————————————————— |
| SWC-0141 | Flexible Silver Nitrate Caustic 6" Applicator Stick (100 pack for Clinic) | AMG | AMG Medical 118-395 | Pack of 100 | N | 1 | For approved ETs nurses and clinic bulk orders only. | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Max | Notes | |
|---------------------|---|-------------|-------------------------|------------------------|-----|---------------|---|---|
| SWC-0121 | No-Sting Barrier Wipe Single Use | Skin-Prep | Smith & Nephew 59420600 | Each 1 mL Wipe | Υ | 3 | Helps reduce friction during removal of tapes and films. Can also be used to prepare skin attachment sites for drainage tubes, external catheters, surrounding ostomy sites and adhesive dressings. | NO-STING SKIN-PREP- Transition of the state |
| SWC-0122 | Dermal Curette with 4mm tip | Miltex | Integra 33-54 | Each | Y | | For sharp/surgical debridement and should only be performed by clinicians with this expertise. | 0 |
| SDR-0006 | Cotton Tipped Applicator Plastic Sterile 2/pk | Medline | MDS202095 | 1 Pack of 2 6" Long | Υ | 7 | | |
| SDR-0010 | Sterile Tongue Depressor | AMG Medical | 118-421 | Adult | Υ | 3 | Adult size Individually wrapped, single patient use. | |

WOUND CARE KITS: Prepared kits of commonly used items. Most items in each kit can be ordered individually if an entire kit is not required if an individual CHRIS order code is listed below. Additional items may also be ordered as needed.

| SDR-0024 | Suture Removal Kit | Medline Canada | Medline DYNDCA1002A | Each | Y | | Single patient use; Sterile, includes Wire Forceps, Gauze, Sponge and Scissors | |
|---------------------|---|-------------------|------------------------|----------|-----|---------------|---|--------------------|
| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
| SMK-0042 In | nitial Wound Care Kit: | N/A | Various | Each Kit | Y | 1 | Maximum Quantity: 1 Kit/Patient Admission | |
| SDR-0039 | Non Allergenic Self-Adhesive Dressing Retention Tape 2" | | | 1 | | | Start up kit only. Visiting Nurse to identify and order regular wound care supplies after initial visit No need to order for patients receiving first visit in | |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per Pack) | | | 4 | | | | |
| SDR-0065 | Sponge Synthetic Sterile Non-Woven 4"x4" (2 per | Pack) | | 4 | | | ambulatory nursing clinics. Clinics are expected to | NO IMAGE AVAILABLE |
| SDR-0006 | Cotton Tipped Applicator Plastic Sterile 2/pk | | | 1 | | | have on hand and replenish stock as needed. | NO IMAGE AVAILABLE |
| SIV-0084 | Syringe Luer Lok 30 mL | | | 1 | | | | |
| SDR-0062 | Wound irrigation tip for syringe | | | 1 | | | | |
| SWC-0107 | Absorbent Dressing 4"x6" | | | 2 | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|-----------------------------|------------|----------------------------|----------|-----|---------------|---|-----|
| SWC-0123 Minor Lace | eration Tray Kit - Sterile: | Busse | Cardinal Health 312-751 | Each Kit | Υ | 1 | Maximum Quantity: 1 Kit/Week | |
| Plastic Tray | | | | 1 | | | Disposable single use For sharp/surgical debridement and should only be | |
| Gauze Sponges | 3 | | | 5 | | | performed by clinicians with this expertise. | ap1 |
| Waterproof Ster | rile Field | | | 1 | | | | |
| Fenestrated Wa | aterproof Drape | | | 1 | | | | |
| Absorbent Towe | els | | | 2 | | | | |
| Tissue Forceps | | | | 1 | | | | 60 |
| Needle Holder | | | | 1 | | | | |
| Scissors | | | | 1 | | | | |

NEGATIVE PRESSURE WOUND THERAPY (NPWT)

| CHRIS | Description | Brand Name | Manufacturer | Size | APO | Weekly | Notes |
|------------|-------------|------------|--------------|------|-----|--------|-------|
| Order Code | | | Code | | | Max | |

Negative Pressure Wound Therapy (NPWT): A therapeutic technique indicated for the application of continual or intermittent negative pressure wound therapy to an acute or chronic wound. The device may promote wound healing by the removal of fluids, including wound exudates, irrigation fluids, body fluids and infections materials. The therapy involves the controlled application of sub-atmospheric pressure to the local wound environment, using a sealed wound dressing connected to a vacuum pump. Equipment and supplies are time-limited as per approved policy.

DISPOSABLE SINGLE-USE NPWT: The primary NPWT devices are to be single-use disposable items when clinically appropriate

PICO Negative Pressure Dressing Kits

| SNP-0069 | PICO-7 Negative Pressure Dressing Kit 4"x8" | PICO 7 | Smith & Nephew 66802002 | 4"x8" 10x20cm | Υ | 1 | 2 Dressings in each kit |
|----------|--|--------|----------------------------|--------------------|---|---|-------------------------|
| SNP-0083 | PICO-7 Negative Pressure Dressing Kit 4"x12" | PICO 7 | Smith & Nephew 66802003 | 4"x12" 10x30 cm | Υ | 1 | 2 Dressings in each kit |
| SNP-0084 | PICO-7 Negative Pressure Dressing Kit 6"x8" | PICO 7 | Smith & Nephew 66802006 | 6"x8" 15x20 cm | Υ | 1 | 2 Dressings in each kit |
| SNP-0085 | PICO-7 Negative Pressure Dressing kit 6"x6" | PICO 7 | Smith & Nephew 66802005 | 6"x6" 15x15 cm | Υ | 1 | 2 Dressings in each kit |



VAC Via Therapy System: Delivers continuous negative pressure (-125mmHg or -75mmHg) for moderate to low exudating wounds, grafts and flaps, (<80ml/day), small to medium-sized wounds (up to 250cm³) equal to or smaller than the size of the medium 3M™ V.A.C.® Via Granufoam™ Spiral Dressing.

| | , ., ., ., ., ., . | VAC Via | 3M Canada VIAKIT077D01 | 1 Kit | Y | | Includes Therapy Unit plus 3 Granufoam Medium Spiral Dressings and 2 Canisters |
|----------|---|---------|---------------------------|-----------------|---|---|---|
| SNP-0087 | VAC Granufoam Medium Spiral Dressing Kit | VAC Via | 3M Canada VIAGFM05.S | 14.5x17x1.75 cm | Y | | Includes VAC Granufoam Medium Spiral Dressing, Drapes, No-sting Barrier Wipe and Wound Ruler |
| SNP-0088 | VAC Via 250 mL Canister | VAC Via | 3M Canada VIACAN05 | 250 mL | Y | 1 | 250mL VAC Via Canister with isolyzer |



See ActiVac Dressing Kits compatible with VAC Via: SNP-0017, SNP-0018, SNP-0019, SNP-0020, and SNP-0022

REUSABLE FOAM NPWT: The ActiVAC NPWT system is the secondary treatment methodology for home and community care. ActiVAC is only to be used if the options for disposable single-use NPWT are not appropriate as authorized by the designated Nurses Specialized in Wound, Ostomy and Continence (NSWOC) in alignment with approved policy.

| ETM-1014 | Reusable Foam NPWT Owned Unit | ActiVAC | 3M Canada 340013 | Each | | 12 weeks | Continuous or Intermittent Modes, Intensity: Low, Medium, High Pressure Settings: -25 to -200 mmHg (-3.3 to -26.6 kPa) Rental only allowed if owned unit is not | |
|----------|--------------------------------|---------|---------------------|------|---|----------|---|--|
| ETM-1010 | Reusable Foam NPWT Rental Unit | ActiVAC | 3M Canada 340013 | Each | Υ | patient | available. Rental unit will automatically be substituted with Owned unit when available as appropriate. | |



| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Max | Notes | |
|---------------------|---|-------------------------|--------------------------|--------------------|----------|----------------|--|--------------------------|
| ActiVAC D canister. | ressing Kits: Includes the dressings, | adhesive drapes | and SensaTRAC Pad. | Excess fluid and p | ootentia | ally infection | ous materials are removed from the wound through the Sen | sa TRAC pad and into the |
| SNP-0008 | VAC GranuFoam Standard Dressing Kit Small | VAC GranuFoam | 3M Canada M8275051/10 | 10x7.5x3.3 cm | Υ | 3 | ActiVAC only: Includes 1 small black foam dressing, 1 adhesive drape,1 SensaTRAC Pad | |
| SNP-0009 | VAC GranuFoam Standard Dressing Kit Medium | VAC GranuFoam | 3M Canada M8275052/10 | 18x12.5x 3.3 cm | Υ | 3 | ActiVAC only: Includes 1 medium black foam dressing, 2 adhesive drapes,1 SensaTRAC Pad | |
| SNP-0010 | VAC GranuFoam Standard Dressing Kit Large | VAC GranuFoam | 3M Canada M8275053/10 | 26x15x3.3 cm | Υ | 3 | ActiVAC only: Includes 1 large black foam dressing, 2 adhesive drapes,1 SensaTRAC Pad | |
| SNP-0066 | VAC GranuFoam Standard Dressing Kit XL | VAC GranuFoam | 3M Canada M8275065/5 | 60x30x1.5 cm | Υ | 3 | ActiVAC only: Includes 1 large black foam dressing, 5 adhesive drapes,1 SensaTRAC Pad. SPECIAL ORDER ITEM - May take up to 5 days for stock. | ()° |
| SNP-0017 | VAC GranuFoam Silver Dressing Kit Small | VAC Silver GranuFoam | 3M Canada M8275098/10 | 10x7.5x3.3 cm | Υ | 3 | ActiVAC and VAC Via: Includes 1 small silver foam dressing, 1 adhesive drape,1 SensaTRAC Pad When using silver dressing kits for VAC Via do not discard the VAC VIA Trac Pad as ActiVAC TRAC Pads and Via TRAC pads are not interchangeable. | |
| SNP-0018 | VAC Granufoam Silver Dressing Kit Medium | VAC Silver Granufoam | 3M Canada M8275096/10 | 18x12.5x3.3 cm | Υ | 3 | ActiVAC and VAC Via: Includes 1 medium silver foam dressing, 2 adhesive drapes,1 SensaTRAC Pad When using silver dressing kits for VAC Via do not discard the VAC VIA Trac Pad as ActiVAC TRAC Pads and Via TRAC pads are not interchangeable. | I or or or or or or or |
| SNP-0019 | VAC Granufoam Silver Dressing Kit Large | VAC Silver Granufoam | 3M Canada M8275099/10 | 26x15x3.3 cm | Υ | 3 | ActiVAC and VAC Via: Includes 1 large silver foam dressing, 2 adhesive drapes,1 SensaTRAC Pad When using silver dressing kits for VAC Via do not discard the VAC VIA Trac Pad as ActiVAC TRAC Pads and Via TRAC pads are not interchangeable. | |
| SNP-0058 | VAC WhiteFoam Small (Foam Only) | VAC Whitefoam | 3M Canada M6275033/10 | 10 x 7.5 x 1 cm | Υ | 3 | ActiVAC and VAC Via 1 polyvinyl alcohol (PVA) foam dressing (10 x 7.5 x 1 cm) | |
| SNP-0059 | VAC WhiteFoam Large (Foam Only) | VAC Whitefoam | 3M Canada M6275034/10 | 10 x 15 x 1 cm | Υ | 3 | ActiVAC and VAC Via 1 polyvinyl alcohol (PVA) foam dressing (10 x 15 x 1 cm) | |
| SNP-0020 | VAC WhiteFoam Dressing Kit Small | VAC Whitefoam | 3M Canada M8275068/10 | 10x7.5x1 cm | Υ | 3 | ActiVAC only: Less adherent material helps promote graft take. Includes 1 small PVA foam dressing, 1 adhesive drape, 1 SensaT.R.A.C. Pad | |
| SNP-0022 | VAC WhiteFoam Dressing Kit Large | VAC Whitefoam | 3M Canada M8275067/10 | 10x15x1 cm | Υ | 3 | ActiVAC only: Less adherent material helps promote graft take. Includes 1 large PVA foam dressing,1 adhesive drape, 1 SensaT.R.A.C. Pad | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---------------------------------------|-----------------|--------------------------|---------------------|----------|---------------|--|----------|
| SNP-0046 | VAC Simplace Dressing Kit Small | VAC Simplace | 3M Canada M8275046/5 | Each | Υ | 3 | ActiVAC only: Includes 1 small Spiral-cut VAC GranuFoam Dressing, 2 adhesive drapes, 1 SensaT.R.A.C. Pad | |
| SNP-0048 | IVAC: Simplace Dressing Kit Medium | VAC Simplace | 3M Canada M8275045/5 | Each | Υ | 3 | ActiVAC only: Includes 1 medium Spiral-cut VAC GranuFoam Dressing, 2 adhesive drapes, 1 SensaT.R.A.C. Pad | |
| ActiVAC A | ccessories: Additional dsposable comp | onents to optim | ise individual treatment | to or when treating | g multip | ole wounds | S. | |
| SNP-0001 | VAC Adhesive Drape | VAC | 3M Canada M6275009/10 | 30.5 x 26 cm | Y | 1 | Included in the Standard Dressing Kits below. This is an adhesive layer spread across the foam dressing, ensuring a closed, moist wound healing environment. | HCCI MAY |
| SNP-0007 | T.R.A.C. Y Connector | VAC | 3M Canada M6275066/10 | Each | Υ | 3 | Y-connector to connect two or more dressings | |
| SNP-0052 | 300 mL ActiVAC Canister with Gel | ActiVAC | 3M Canada M8275058/10 | Each | Υ | 1 | 1 canister with gel, tubing, clamp and connector; May be authorized for more than 1/week if the canister is full. | |
| SNP-0060 | VAC SensaT.R.A.C. Pad | VAC | 3M Canada M8275057/10 | Each | Y | I I | Includes 1 SensaT.R.A.C. Pad with tubing, clamp and connector | |
| SNP-0061 | T.R.A.C Tubing Cap | VAC | 3M Canada M6275069/10 | Each | Υ | 1 | Cap to secure end of tubing when not in use | - |
| SNP-0067 | VAC GranuFoam Bridge Dressing Kit | VAC Bridge | 3M Canada M8275042/5 | Each | Υ | 3 | Each kit includes black foam, bridge dressing; sensa TRAC pad and adhesive drape. SPECIAL ORDER ITEM | |
| SNP-0068 | VAC GranuFoam Bridge XG Dressing Kit | VAC Bridge | 3M Canada M8275044/5 | Each | Υ | 3 | - May take up to 5 days for stock. | |

CLEANSERS & ANTISEPTIC SOLUTIONS

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Max | Notes |
|---------------------|--|--------------------|-----------------------------|-----------------|-----|------------------|---|
| Cleansers | & Antiseptic Solutions | | | | | | |
| SSO-0001 | Chlorhexidine 2% Alcohol 70% 500 mL Bottle | SoluPrep | 3M 103.15 | 500 mL | Υ | ' | Reclosable bottle. Document date opened on bottle and discard unused product after 30 days |
| SSO-0002 | Hydrogen Peroxide 3% 450mL Bottle | Swan or equivalent | Various | 450-500 mL | Υ | 1 | Reclosable bottle. Document date opened on bottle and discard unused product after 7 days |
| SSO-0019 | Normal Saline 500 mL bottle | Baxter | Baxter JF7123 | 500 mL | Υ | 1 | When used for for wound care these items are approved only for: 1) Immunocomprimised patient (i.e. BMT, Oncology, Organ |
| SSO-0023 | Normal Saline 118 mL Irrigation Bottle Sterile | Derma | Derma Sciences DUP48501N | 118 mL | Υ | 7 | Transplant patient); 2) Well water (untreated) or |
| SSO-0031 | Sterile Water 500 mL bottle | Baxter | Baxter JF7113 | 500 mL | Υ | 1 | 3) Environmental considerations (i.e. concerns of hygiene) |
| SSO-0069 | Chlorhexidine 2% Alcohol 70% Wipe | SoluPrep | 3M 101.06 | Each | Υ | 28 | Single use antiseptic skin preparation wipe for swabbing sites and ports for vascular access catheters. Single use |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | SoluPrep | 3M 102.03 | Small 1.6 mL | Υ | 8 | Antiseptic skin preparation swabsticks for maintenance of indwelling vascular access catheters. Single use |
| SSO-0052 | Alcohol 70% Wipe | -Cardinal | Loris | 30 x 65 mm | Υ | 28 | Single use antiseptic topical skin cleanser prior to an injection Medium-sized wipe 70% Isopropyl Alcohol; Single use |
| SSO-0072 | Alcohol 70% Wipe - Box of 200 | Carumai | LP103-03 | 30 x 03 IIIII | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SSO-0054 | Povidone Iodine 10% Wipe | Cardinal | Loris LP109-01 | 30 x 65 mm | Υ | 3 | Single use 10% PVP-I antiseptic/germicide prep used for wound clensing or minor invasive procedures, scrapes, cuts or burns. Non- |
| SSO-0067 | Povidone Iodine 10% Swabstick (Single) | Cardinal | Loris LP109-05 | 13 x 25 mm | Υ | 8 | irritating and virtually non-stinging for better patient comfort. |
| SSO-0004 | Povidone Iodine 10% 115mL Bottle | Loris | Loris LP109-09 | 115 mL | Υ | 2 | Review manufacturer's product information regarding contraindications. |
| SSO-0011 | Hygeol 1:20 500 mL | N/A | N/A | 500 mL | Y | 7 per 2 wks | Maximum 14/patient |
| SSO-0073 | Wound Irrigation Gel Antimicrobial 30 mL Bottle | Prontosan | B. Braun BB400510 | 30 mL | Υ | 1 per 28 days | For cleansing, moistening and decolonization of acute and chronic skin wounds, 1st and 2nd degree burns. (NERDS: Non-healing, |
| | Wound Irrigation Solution Antimicrobial 350 mL Bottle | Prontosan | B, Braun BB400431 | 350 mL | Υ | 1 per | Exudate, Red/bleeding, Debris, Smell). This product is not to be used to replace potable water cleansing |
| Miscellane | ous | | | | | - | _ |
| SO1-0083 | Oral Swab sticks Latex-free Plain Untreated (Box 250) | Toothette | Sage Products 5602UT | Box of 250 | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SOT-0084 | Oral Swab sticks Latex-free Plain Untreated (10 Pack) | Toothette | Sage Products 6072 | Pack of 10 | Υ | 2 | Palliative patients only at home/in community |

I.V. SUPPLIES

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|--------------|----------------------|--------------|-----|----------------------|---|
| NEEDLES 8 | SYRINGES | | | | | | |
| SIV-0079 | Syringe Luer Lok 1mL | Luarilak | BD | 41 | Υ | 7 | Luer Lok tip - Graduation 1/100 mL |
| SIV-0238 | Syringe Luer Lok 1mL - Box of 100 | -Luer-Lok | 309628 | 1 mL | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0080 | Syringe Luer Lok 3 mL | Luor Lok | BD | 0 ml | Υ | 7 | Luer Lok tip - 0.2 mL graduations |
| SIV-0239 | Syringe Luer Lok 3 mL - Box of 200 | -Luer-Lok | 309657 | 3 mL | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0081 | Syringe Luer Lok 5 mL | Luer-Lok | BD 309646 | 5 mL | Υ | 7 | Luer Lok tip - Graduation 1/5 mL |
| SIV-0082 | Syringe Luer Lock 10 mL | Luer-Lok | BD 302995 | 10 mL | Υ | 7 | Luer Lok tip - Graduation 1/5 mL |
| SIV-0084 | Syringe Luer Lok 30 mL | Luer-Lok | BD 302832 | 30 mL | Υ | See Notes | Luer Lok tip - Graduation: 1 mL For IV: 7 / Patient / Week For Wound Care: 1 / Patient / Week In-Home or 1 / Visit in Clinic Also order Wound Irrigation Tip for syringe for each dressing change (SDR-0062) |
| SIV-0085 | Syringe Luer Lok 50 mL | Luer-Lok | BD 309653 | 50 mL | Υ | 7 | Luer Lok tip - Graduation: 1 mL , 2 oz in 1/4 oz |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | DD. | BD | 18g 1.5" | Υ | 3 | |
| SIV-0240 | Blunt Fill Needle 18g 1.5" - Box of 100 | BD | 305180 | 1.2 x 40 mm | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0173 | Needle with Syringe 3 mL 23g 1" | Eclipse | BD 305782 | 3mL 23g 1" | Υ | 7 | Eclipse Needle with Luer-Lok |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | Eclipse | BD 305780 | 1mL 25g 5/8" | Υ | 7 | Eclipse Needle with Luer-Lok |
| SIV-0177 | Needle with Syringe 1mL 27g 0.5" | Eclipse | BD | 1ml 27g 0 5" | Υ | 7 | Eclipse Needle with Luer-Lok |
| SIV-0241 | Needle with Syringe 1mL 27g 0.5" - Box of 50 | T L C II PSE | 305789 | 1mL 27g 0.5" | N | 3 Boxes per Month | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0181 | Needle with Syringe 3mL 25g 5/8" | Eclipse | BD 305781 | 3mL 25g 5/8" | Υ | 7 | Eclipse Needle with Luer-Lok |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|--------------------------|------------------------------|----------------------|-------|-------------------|---|
| SIV-0192 | Gripper Plus 22g x 3/4" Non Y-site | Gripper Plus | Smiths Medical 21-2766-24 | 22 gauge x 3/4" | Υ | 1 | Non-coring needle without Y-site Used to access implanted vascular access devices. |
| SIV-0193 | Gripper Plus 19g x 1.25" Non Y-site | Gripper Plus | Smiths Medical 21-2765-24 | 19 gauge x 1.25" | Υ | | Non-coring needle without Y-site Used to access implanted vascular access devices. |
| SIV-0211 | Blunt Fill Needle with Filter 18g x 1.5" | ·BD | BD | 18g 1.5" | Υ | | 5 micron filter with blunt fill tip. Only to be used in Symptom Relief Kits for use when withdrawing fluid/medication from an ampoule |
| SIV-0242 | Blunt Fill Needle with Filter 18g x 1.5" - Box of 100 | טם | 305211 | 1.2 x 40 mm | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| PREFILLED | SYRINGES: For the sole purpose of flushing | Vascular Acce | ss Devices as per | nursing guideline | s and | manufactur | er's recommendations. Home and Community Care Support |
| | rth East shall not fund or provide any other co | ncentrations of | • | 0 0 | | | , |
| SIV-0196 | Heparin Lock Syringe 10 USP units/mL - 5 mL Fill in 10 mL Diameter Syringe | PosiFlush | BD 306414 | 10 USP/mL 5 mL/mL | Υ | 3 | |
| SIV-0199 | Heparin Lock Syringe 10 USP units/mL - 3 mL Fill in 10 mL Diameter Syringe | PosiFlush | BD 306413 | 10 USP/mL 3 mL | Υ | 3 | |
| SIV-0197 | Heparin Lock Syringe 100 USP units/mL - 5 mL Fill in 10 mL Diameter Syringe | PosiFlush | BD 306424 | 100 USP/mL 5 mL | Υ | 3 | |
| SIV-0200 | Heparin Lock Syringe 100 USP units/mL - 3 mL Fill in 10 mL Diameter Syringe | PosiFlush | BD 306423 | 100 USP/mL 3 mL | Υ | 3 | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | PosiFlush | BD 306592 or equivalent | 10 mL | Υ | 14 | |
| I.V. SOLUTI | ONS | | | | | | |
| SSO-0024 | Sodium Chloride 0.9% INJ USP 100 mL Bag | Viaflex | Baxter JB1309M | 100 mL | Υ | 7 | |
| SSO-0025 | Sodium Chloride 0.9% INJ USP 250 mL Bag | Viaflex | Baxter JB1322P | 250 mL | Υ | 7 | |
| SSO-0026 | Sodium Chloride 0.9% INJ USP 500 mL Bag | Viaflex | Baxter JB1323 | 500 mL | Υ | 7 | |
| SSO-0027 | Sodium Chloride 0.9% INJ USP 1000 mL Bag | Viaflex | Baxter JB1324 | 1000 mL | Υ | 7 | |
| SSO-0029 | Sodium Chloride 0.9% 10 mL Vial Preservative-free | Hospira or equivalent | 954354 | 10 mL | Υ | | Indicated only for diluting or dissolving medications for use in a nebulizer |
| SSO-0032 | Sterile Water 10 mL Vial | Hospira or equivalent | 954469 | 10 mL | Υ | , | Indicated only for diluting or dissolving medications for use in a nebulizer |
| SSO-0070 | Dextrose 5% & Sodium Chloride 0.45% INJ USP 1000 mL Bag | Viaflex | Baxter JB1074 | 1000 mL | Υ | 7 | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|--------------------------|------------------------------|---|-----|-------------------|--|
| SSO-0034 | Dextrose 3.3% & Sodium Chloride 0.3% INJ USP 1000 mL Bag | Viaflex | Baxter JB1034 | 1000 mL | Υ | 7 | |
| SSO-0041 | Dextrose 5% INJ USP 1000 mL Bag | Viaflex | Baxter JB0064 | 1000 mL | Υ | 7 | |
| SSO-0071 | 10% Dextrose INJ USP 1000 mL Bag | Viaflex | Baxter JB0164 | 1000 mL | Υ | 7 | |
| SSO-0046 | Lactated Ringer's INJ USP 1000 mL Bag | Viaflex | Baxter JB2324 | 1000 mL | Υ | 7 | |
| I.V. CATHET | TERS & ACCESS | | | | | | |
| SIV-0150 | IV Catheter 22g Blue | Nexiva | BD 383512 | 22 gauge1.00" 0.9 mm x 25 mm | Υ | 3 | Closed IV catheter system-Single port |
| SIV-0151 | IV Catheter 24g Yellow | Nexiva | BD 383511 | 24 gauge 0.75" 0.7 mm x 19 mm | Υ | 3 | Closed IV catheter system-Single port |
| SIV-0228 | IV Catheter 18g Green | Nexiva | BD 383519 | 18 gauge 1.25" 1.1 mm x 32 mm | N | | Closed IV catheter system-Single port For MAiD Patients as approved by Care Coordinator |
| SIV-0229 | IV Catheter 20g Pink | Nexiva | BD 383517 | 20 gauge 1.25" 1.3 mm x 32 mm | N | // | Closed IV catheter system-Single port For MAiD Patients as approved by Care Coordinator |
| SIV-0187 | IV Catheter Set - 22g | Saf-T-Intima | BD 383329 | 22 gauge 0.75" 0.9 mm x 19 mm | Υ | 3 | Safety system with wings, extension tubing and Y adapter |
| SIV-0203 | IV Catheter Set - 24g | Saf-T-Intima | BD | 24 gauge 0.75" | Υ | 3 | Safety system with wings, extension tubing and Y adapter |
| SIV-0243 | IV Catheter Set - 24g - Box of 25 | oar i illillia | 383319 | 0.9 mm x 19 mm | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| I.V. GRAVIT | Y SETS AND EXTENSION SETS | | | | | | |
| SIV-0225 | Continuous Flow IV Solution Set | Clearlink Continu-Flo | Baxter 2C8519 | 10 drops / mL 112" (2.8 m) | Υ | / | 2 Luer Activated Valves, Male Luer Lock Adapter, Approximate Volume 10 drops per mL |
| SIV-0226 | Secondary Medication Set | Clearlink | Baxter JC7453 | 10 drops / mL 37" (94 cm) | Υ | 7 | With Male Luer lock adapter and hanger |
| SIV-0009 | Needle-Free Connector Valve | MaxZero | Cardinal Health MZ1000-07 | Length: 3 cm Priming Volume 0.19 mL | Y | 7 | Clear needleless connector |
| SIV-0012 | Catheter Extension Set | MaxZero | Cardinal Health MZ5301 | Length: 18 cm Priming Volume 0.3 mL | Υ | 3 | Extension set with clear needle-free connector valve |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|------------------------------|------------|---------------------------|-------------------------|-----|-------------------|---|
| SIV-0224 | Y Connector IV Extension Set | MaxZero | Cardinal Health MZ5307 | 7" / 8 cm PV: 0.8 mL | Υ | 1 | Bifuse pressure rated extension set, 2 bonded MaxZero connectors, 2 slide clamps, spin male luer lock |

I.V. SECUREMENT: Securement devices and dressings used to maintain catheter securement, and protect peripheral and central catheter insertion sites and wounds. The primary product is the Tegaderm IV Advanced and IV3000 is only to be authorized for patient-specific cases if Tegaderm IV dressings are clinically not appropriate due to sensitivity.

| SDR-0094 | Tegaderm IV Advanced Securement Dressing | Tegaderm IV | 3M 1683 | 2.5 "x 2.75" 6.5 cm x 7 cm | Υ | 3 | For peripheral IVs |
|----------|--|-------------------------|------------------------|----------------------------------|---|---|---|
| SIV-0244 | Tegaderm IV Advanced Securement Dressing - Box of 100 | regademirv | | | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0231 | PICC/CVC Securement Device and Dressing Kit (Regular) | Tegaderm IV Advanced | 3M 1837-2100 | 3 .5" x 4.5" 8.5 cm x 11.5 cm | V | 1 | To secure Peripherally Inserted Central Catheters (PICC's) and short term Central Venous Catheters (CVC's) Designed to be worn for up to 7 days Removes cleanly, without causing patients undue pain or distress Accommodates the majority of PICC and CVC catheters up to and including 12 French |
| SIV-0232 | PICC/CVC Securement Device and Dressing Kit (Large) | Tegaderm IV Advanced | 3M 1839-2100 | 4" x 6" 10 cm x 15.5 cm | ı | | |
| SIV-0236 | IV3000 Securement Dressing | IV 3000 1Hand | Smith & Nephew 4008 | 4" x 4-3/4" 10cm x 12cm | Υ | | Only to be ordered if Tegaderm IV dressings are not appropriate. Dressing only for Central/PICC/midline catheters. |
| SIV-0134 | StatLock IV Catheter Securement Device | StatLock PICC Plus | Bard PIC0220 | Each | Y | 1 | Secures catheter in place when using the IV3000 dressing when required to secure PICC line. Accommodates a wide variety of PICC Catheters. Tricot Anchor Pad, Sliding Post, Sterile. |

NOTE: Additional items are also available in the TRACHEOTOMY / DRAIN SPONGES in the Dressings, Bandages, Tape, Compression & Wraps section

INFUSION PUMPS & SUPPLIES

| ETM-1005 | IV pole - portable | Rental equipment | |
|----------|---|------------------|---|
| ETM-1009 | Ambulatory infusion pump (includes pump, power supply, small carrying pouch and batteries) | Rental equipment | All service areas |
| ETM-1012 | Large carry case - for ambulatory infusion pump (if required)(For infusions requiring 250 mL or 500 mL bags) | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| ETM-1016 | Remote Dose Extension Cord for ETM-1009 IV Pump | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| ETM-0028 | Sapphire IV Pump Mini-Cradle Mount for Pole | Rental equipment | Cochrane sub-region/Timmins Branch area (Bayshore Specialty Rx) |
| ETM-0030 | Sapphire IV Pump Wearable Pouch 500 ml | Rental equipment | Cochrane sub-region/Timmins Branch area (Bayshore Specialty Rx) |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|---|------------|------------------------------|----------------|-----|-------------------|---|
| ETM-0034 | Sapphire IV Pump Bolus Handle | | Rental e | equipment | | | Cochrane sub-region/Timmins Branch area (Bayshore Specialty Rx) |
| SOT-0026 | Replacement batteries for CADD IV pump | | | 4- AA's | Y | | Single use batteries for CADD infusion pump- Solis Not used in Timmins Branch/Cochrane sub-region |
| | | Bodyguard | CME BBA120-161XPS | 80" | Υ | 7 | Timming Propoly/Cookrang out region (Poughers Specialty Dy) |
| SIV-0042 | Ambulatory Infusion Pump Administration Set | Sapphire | AP-407 12000-000-0008 | 112" 286 cm | Υ | 7 | Timmins Branch/Cochrane sub-region (Bayshore Specialty Rx) |
| | | CADD | Smiths Medical 21-7022-24 | Each | Υ | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0043 | Ambulatory Infusion Pump High Volume Administration Set | CADD | Smiths Medical 21-7345-24 | Each | Υ | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| | | Bodyguard | CME A120-112XSXFK | 110" | Υ | 7 | Timmins Branch/Cochrane sub-region (Bayshore Specialty Rx) |
| SIV-0044 | Ambulatory Infusion Pump Administration Set Filtered | Sapphire | AP 206 12000-000-0005 | 116" 295 cm | Υ | 7 | |
| | | CADD | Smiths Medical 21-7346-24 | Each | Υ | / | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0045 | Ambulatory Infusion Pump Extension Tubing Set | ICU | ICU Medical B1327 | 20 inches | Y | 3 | Timmins Branch/Cochrane sub-region (Bayshore Specialty Rx) |
| 317-0045 | Ambulatory initiasion Fump Extension Tubing Set | CADD | Smiths Medical 21-7045-24 | 30 Inches | Ţ | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0046 | Ambulatory Infusion Pump Extension Tubing 45 inches | CADD | Smiths Medical 21-7046-24 | 45 inches | Υ | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0047 | Ambulatory Infusion Pump Medication Cassette 50ml | | | 50 mL | Υ | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0048 | Ambulatory Infusion Pump Medication Cassette 100ml | | | 100 mL | Υ | 7 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0202 | IV Medication Refill Solution Bag | | | Each | Υ | 7 | Solution bag - not a cassette. For IV Med Refill orders only |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | | | | | |
|---------------------|--|------------------------|------------------------|-------------------------|--------|-------------------|--|--|--|--|--|--|
| physician, el | INGLE USE ELASTOMERIC DISPOSABLE PUMP & SUPPLIES: The primary IV pump is a re-usable ambulatory infusion pump. If a single-use pump is requested by the patient or hysician, elastomeric/single-use pumps shall only be provided if medications are funded under the Ontario Drug Benefit (ODB). Other costs related to elastomeric pumps that are not covered y ODB are the sole responsibility of the patient. See also <i>Infusion Therapy Policies</i> | | | | | | | | | | | |
| SIV-0233 | Elastomeric Single Use Disposable IV Pump | EasyPump or equivalent | B.Braun | Each | Ν | 1 | SPECIAL ORDER ITEM: Care Coordinator to confirm ODB coverage and timing for availability with the approved IV pharmacy prior ordering pump. The pump type/size is selected by pharmacy per medication type. | | | | | |
| SIV-0234 | Elastomeric IV Pump Catheter Extension Set | Max Plus Clear | MP9014C | 14"/36 cm PV 0.75 ml | Υ | 3 | For use with SIV-0233 Elastomeric pump only. | | | | | |
| SIV-0235 | Elastomeric IV Pump Belt Carrying Bag | EasyPump or equivalent | B.Braun | Each | N | 1 / patient | For use with SIV-0233 Elastomeric pumps only. | | | | | |
| BLOOD CO | LLECTION SUPPLIES: For blood collection for | or patients eligil | ole for in-home Lal | b Services as initi | ated b | y Care Coc | ordinator only as per policy | | | | | |
| SIV-0137 | Stretch Tourniquet | Vacutainer | BD 367203 | 18" x 1" | Υ | 1 | | | | | | |
| SIV-0145 | Blood collection infusion set with luer adapter 23g | Safety-Lok | BD 367283 | 23 g x 0.75" Needle | Υ | 3 | 12" Tubing Manually activated safety shield to fully cover needle. | | | | | |
| SIV-0204 | Blood collection infusion set with luer adapter 21g | Safety-Lok | BD 367281 | 21 g x 0.75" Needle | Υ | 3 | 12" Tubing Manually activated safety shield to fully cover needle | | | | | |
| SIV-0212 | Blood Tube holder Single Use Nonstackable | Vacutainer | BD 364815 | Each | Υ | 1 | Single patient use only - Discard the entire device. Needle is not part of the tube holder. One Blood Tube Holder can be used for all the blood work per patient | | | | | |
| SIV-0205 | Venous Blood Collection Tube - Lavender | Vacutainer | BD 367841 | 13 x 75 mm x 2.0 mL | Υ | 3 | | | | | | |
| SIV-0206 | Venous Blood Collection Tube - Yellow | Vacutainer | BD 367986 | 13 x 100 mm x 5.0 mL | Υ | 3 | | | | | | |
| SIV-0207 | Venous Blood Collection Tube - Blue | Vacutainer | BD 363083 | 13 x 75 mm x 2.7 mL | Υ | 3 | | | | | | |
| SIV-0208 | Venous Blood Collection Tube - Green | Vacutainer | BD 367962 | 13 x 100 mm x 4.5 mL | Υ | 3 | | | | | | |
| SIV-0209 | Venous Blood Collection Tube - Red | Vacutainer | BD 367815 | 13 x 100 mm x 6.0 mL | Υ | 3 | | | | | | |
| SIV-0210 | Lab Specimen Bag | Fisher Scientific | Fisherbranded 22310004 | 6" x 10" | Υ | 1 | Transparent Blue, Non-Sterile Single Use Disposable With separate outside pocket holds the requisition slip | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|---------------------|--------------------------------------|------------------------------|-----|-------------------|--|
| OTHER I.V. | SUPPLIES & ACCESSORIES | | | | | | |
| SIV-0132 | Sterile non-fenestrated drape 18x26" | Busse or equivalent | Cardinal Health 696 or equivalent | 18 x 26" | Υ | 2 | Multi-layered, absorbent, all-purpose sterile field |
| SIV-0213 | Vacuum Drainage Bottle 1000mL with 120cm Tubing Set | PreVac MedRx | CHS S55148 | 1000 mL | Y | 3 | Sterile bottle includes attached drain line with male luer connector and roller camp (drain at 12 PSI). Not to be used for body fluid samples for diagnostic testing |
| SIV-0214 | Drainage Procedure Kit (Excludes Bottle) | MedRx | CHS 2976C | 1 Kit | Υ | 3 | Order SIV-0213 drainage bottle separately. Includes dressings, tubing, alcohol wipes, vented luer cap, tube clamp, and chemo-rated Nitrile gloves. Sterile. |
| SIV-0215 | Sterile Vented Luer Cap | LSL Health Care | CHS 797A | Each | Υ | 3 | Not required if ordering Kit 2976C (SIV-0214) For pts with luer lock connections |
| SIV-0245 | Drainage Bag 2000mL with 120cm Tubing Set | MedRx | CHS M04200 | 2000 mL | Υ | 7 | Includes 120 cm detachable drain line with rotating male luer fitting; roller clamp; built in rotating hanger; integrated drip chamber and anti-reflux valve |
| SIV-0217 | Safety Drain Line Adapter | MEDRx | CHS SA5200 | Each | Υ | 3 | For use ONLY if patient has Pleurex™ catheter |
| SIV-0227 | Universal Catheter Cap | Medionics | CHS S55300 | Each | Υ | 3 | For use ONLY if patient has Pleurex™ catheter |
| SIV-0131 | Sharps Container Phlebotomy | SharpSafety | Cardinal Health Z8900SA | 6.25x4.25x4.25" 1 Litre | Υ | 1 | Accommodates a variety of Sharps; Temporary closure for security between uses; Lids lock for safe transport and disposal. |
| SIV-0218 | Sharps Container Chemotherapy | ChemoSafety | Cardinal Health 8982CRH | 10"x7.25"x10.5" 7.6 Litre | Υ | 1 | For the disposal of chemotherapy or other cytotoxic waste or for MAiD patients as approved by the Care Coordinator. Temporary closure and locking lid for safe transport and disposal. |
| SIV-0219 | Spill Kit for Cytotoxic Drugs | Safetec | Safetec 14399 | Each | Υ | 1 | Protects and cleans up spills of chemo or other hazardoudous medications. All inclusive kits contains all needed supplies and instructions. |

| CHRIS Order Code | | Brand Name | Manufacturer Code | Size | АРО | Maximum | Notes | | | |
|---------------------|---|------------------|------------------------------|---------------|--|---------|---|--|--|--|
| | epared kits of commonly used items. Most iter llso be ordered as needed. | ns in each kit c | re kit is | s not require | ed and the individual CHRIS order code is listed below. Additional | | | | | |
| | ICC Dressing Kit: | 1 | Maximum Quantity: 1 Kit/Week | | | | | | | |
| SDR-0005 | Dressing Tray Kit | | | 1 | | | All PICC/central lines require pump for infusion therapy. A large Ambulatory Infusion Pump carrying bag may be required, if necessary. | | | |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1 m | | | Add: 2 of Ambulatory Infusion Pump Administration Sets (SIV-0042) | | | |
| SIV-0231 | PICC/CVC Securement Device and Dressing Kit (| Regular) | | 1 | | | OR Ambulatory Infusion Pump Administration Sets Filtered for administration of Cephalosporins, Vancomycin, Ceftazidine (SIV-0044) | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | and 2 Alcohol Wipes (SS0-0052) AND 2 Sodium Chloride 0.9% Syringes (SSO-0021) per day of infusion | | | |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask 2 | | | | | | Add: Portacath Gripper item SIV-0192 or SIV-0193 as required. | | | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | | | 1 | | | Add: Prefilled Heparin Lock Syringes as per nursing guidelines (See IV, Needles, Syringes) | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 2 | | | Dressing Tray SDR-0005 may be ordered individually if this entire PICC kit is not required. | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 2 | | | | | | |
| For double l | umen add the following to your order: | | | | | | | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | Υ | | | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 1 | Υ | | | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 2 | Υ | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO Y | Weekly Maximum | Notes |
|--|---|------------|----------------------|------|----------|--|---|
| SMK-0039 Subcutaneous Administration of Pain Management Medications via Infusion Pump: Initial Start Up Kit | | | | | | | Maximum Quantity: 1 Kit/Week |
| SWC-0072 | Transparent Film Dressing 10 x 12cm | | | 2 | 1 | | Add: Extension tubing (SIV-0045 or SIV-0046) as needed if being used with medication cassette. |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | 1 | | |
| SIV-0203 | IV Catheter Set - 24g | | | 2 | 1 | | |
| SIV-0009 | Needle-Free Connector Valve | 1 | | | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | 4 | | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 1 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 4 | | | |
| SMK-0040 P | eripheral IV Start Kit: | | | | Υ | 2 | Maximum Quantity: 2 Kits/Week |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | | | 2 | 1 | | For infusion via pump order: Ambulatory Infusion Pump Administration Set (SIV-0042) |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1 m | 1 | of Cephalosporins, Vancomycin, Ceftazidine (SIV-0044) <u>AND</u> | Ambulatory Infusion Pump Administration Set Filtered for administration |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | AND |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | 2 Chlorhexidine 2% Alcohol 70% Wipes (SS0-0069) AND |
| SIV-0137 | Stretch Tourniquet | | | 1 | | | 2 Sodium Chloride 0.9% Syringes (SSO-0021) per day of infusion |
| SIV-0150 | IV Catheter 22g Blue | | 2 | | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | 4 | | | | |
| SSO-0052 | Alcohol 70% Wipe | | 4 | | | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | 2 | | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 2 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | | |
|---------------------|---|------------|----------------------|------|------------------------------|-------------------|---|--|--|
| SMK-0041: 0 | Gravity Peripheral IV Administration Kit | | Υ | 2 | Maximum Quantity: 2 Kit/Week | | | | |
| SIV-0137 | Stretch Tourniquet | | | 1 | - | | Provides 3 days of infusion For increased frequency of medication administration beyond once | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | daily, order additional Alcohol wipes (SSO-0052) and Prefilled Sodium Chloride 0.9% Syringes (SSO-0021) per day of infusion | | |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | | | 2 | | | If additional line required add SIV-0226 Secondary IV set to the order | | |
| SIV-0225 | Continu-Flo solution set | | | 3 | | | Additional catheters may be ordered as needed | | |
| SSO-0024 | Sodium Chloride 0.9% INJ USP 100 mL Bag | | | 3 | | | Additional catheters may be ordered as needed | | |
| SIV-0151 | IV Catheter 24g Yellow | | | 2 | | | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | | | |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1 m | | | | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 1 | | | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 6 | | | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 6 | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|------------|----------------------|------|--|-------------------|--|
| | RELIEF SUPPLY KITS: Geographic area-spe h kit can be ordered individually if an entire kit | | | | RC 95 End-of-Life patients. Prepared kits of commonly used items. | | |
| SIV-0223: Syı | mptom Relief Supply Kit - North Bay | | Υ | 1 | Medication order must be provided to the geographically designated approved Infusion Supplies pharmacy provider with the supplies order. The approved Infusion Supplies provider will deliver medication and | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | supply kit to the patient's home or alternate treatement address specified in the order |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per | Pack) | | 2 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 2 | | | |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | | 5 | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 10 | | | Medication order must be provided to the pharmacy with the supplies order. Includes home delivery within Greater Sudbury only: Designated approved Infusion Supplies pharmacy provider will deliver |
| SIV-0203 | IV Catheter Set - 24g | | | 2 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 20 | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 2 | | | |
| SIV-0220: Syı | mptom Relief Supply Kit - Sudbury | | | | Y | 1 | |
| SIV-0131 | Sharps Container Phlebotomy | | | 1 | | | medication and supply kit to the patient's treatement address as specified in the order. |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | specified in the order. |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 4 | | | |
| SIV-0203 | IV Catheter Set - 24g 4 | | | | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 25 | 1 | | |
| SSO-0052 | Alcohol 70% Wipe | | | 20 |] | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|----------------------|-------------------|-------|
| SIV-0221: Sy | mptom Relief Supply Kit - Timmins | | Υ | 1 | BSRx Item No: 300301 | | |
| SSO-0052 | Alcohol 70% Wipe | | 25 | | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 25 | | | |
| SIV-0080 | Syringe Luer Lok 3 mL | | | 15 | | | |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | | 40 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 4 | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 4 | | | |
| SIV-0131 | Sharps Container Phlebotomy | | 1 | | | | |
| SIV-0009 | Needle-Free Connector Valve | | 4 | | | | |
| SIV-0230: Sy | mptom Relief Supply Kit - Kirkland Lake | | | | Υ | 1 | |
| SSO-0052 | Alcohol 70% Wipe | | | 25 | | | |
| SIV-0211 | Blunt Fill Needle with Filter 18g x 1.5" | | | 7 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 4 | | | |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per | Pack) | | 4 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 40 | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | 20 | - | | | | |
| SIV-0203 | IV Catheter Set - 24g | 4 | | | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | 4 | | | | | |
| SIV-0131 | Sharps Container Phlebotomy | 1 | | | | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|------------|----------------------|------|-----|-------------------|--|
| SIV-0237: Sy | SIV-0237: Symptom Relief Supply Kit - Manitoulin | | | | | 1 | The medication order may be provided by any local pharamcy, but the supplies are provided by the designated approved Infusion Supplies pharmacy provider |
| SSO-0052 | Alcohol 70% Wipe | | | 25 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 25 | | | |
| SIV-0080 | Syringe Luer Lok 3 mL | | | 15 | | | |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | | 40 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 4 | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 4 | | | |
| SIV-0131 | Sharps Container Phlebotomy | | | 1 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 4 | | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 20 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|-----|--|-------|
| SIV-0222 Syı | mptom Relief Supply Kit - Sault Ste Marie | | | Y | 1 | The medication order may be provided by any local pharamcy, but the supplies are provided by the designated approved Infusion Supplies pharmacy provider | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 3 | | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 3 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 3 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 1 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 30 | | | |
| SIV-0117 | Blunt Fill Needle 18g x 1.5:" | | | 30 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 30 | | | |
| SIV-0080 | Syringe Luer Lok 3 mL | | | 5 | | | |
| SIV-0131 | Sharps container Phelbotomy | | | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|---|-------------------|--|
| MAIDKIT: MA | AiD Supplies Kit | | N | 1 | Pre-assembled kit of supplies required for use by the physican with eligible, active patients approved for MAiD as per approved Home and Community Care Support Services North East policy. | | |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | | | 4 | | | When placing an order from Robinson's, add the delivery code #FC- |
| SDR-0004 | Gauze Sterile Woven 4"x4" (2 per Pack) | | | 4 | | | 1033 for LHIN-funded Home Delivery for MAiD to this order. This Finance Charge is in addition to the appropriate Urban or Rural Delivery |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1m | | | Finance Charge and only applicable to this kit. |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 2 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 2 | | | |
| SIV-0012 | Catheter Extension Set 18 cm | | | 2 | | | |
| SIV-0218 | Sharps Container Chemotherapy | | | 1 | | | |
| SIV-0137 | Stretch Tourniquet | | | 1 | | | |
| SIV-0224 | Y Connector IV Extension Set | | | 2 | | | |
| SIV-0228 | IV Catheter 18g Green | | | 2 | | | |
| SIV-0229 | IV Catheter 20g Pink | | | 4 | | | |
| SOT-0096 | Non Sterile Nitrile Gloves Large (Pair) | | | 4 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 8 |] | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 6 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 8 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|----------------|----------------------|------|--|-------------------|--|
| SMK-0045 - I | Remdesivir Kit for Gravity Infusion | | | Qty | APO | Weekly Maximum | Notes |
| SIV-0225 | Continuous Flow IV Solution Set | 3 | | 1 | This supply kit supports up to 3 doses | | |
| SIV-0009 | Needle Free Connector | | | 3 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 18 | | | |
| SSO-0021 | Sodium Chloride 0.9% 10mL prefilled syringe | | | 12 | | | |
| SIV-0210 | Lab Specimen Bag | | | 3 | | | |
| SMK-0046 - I | Remdesivir Kit for Pump Infusion | | | Qty | АРО | Weekly Maximum | Notes |
| SIV-0043 | High Volume Administration Set (or If Bayshore R | x SIV-0042 Pum | np Admin. Set) | 3 | | 1 | This supply kit supports up to 3 doses |
| SIV-0009 | Needle Free Connector | | | 3 | | | If no PIV Insitu then order the kit and add: 1 - SMK-0040 - Peripheral IV Start Kit |
| SSO-0052 | Alcohol 70% Wipe | | | 18 | | | , i |
| SSO-0021 | Sodium Chloride 0.9% 10mL prefilled syringe | | | 15 | | | |
| SIV-0210 | Lab Specimen Bag | | | 3 | | | |
| SMK-0047 - I | Remdesivir Reconstitution Kit 3 day Treatment (| Timmins) | | Qty | АРО | Weekly Maximum | Notes |
| SSO-0032 | Sterile Water 10 mL Vial | | | 8 | | 1 | Timmins/Cochrane Area ONLY This supply kit supports up to 3 doses |
| SSO-0025 | NaCl 0.9% 250 mL Bag | | | 3 | | | For central line add: |
| SIV-0117 | Blunt Fill Needle | 3 | | | 2 of SSO-0006 - CHG 2% ALC 70% swabstick add 2 of SSO-0021 - Normal Saline 10 mL | | |
| SSO-0052 | Alcohol 70% Wipe | | | 18 | | | |
| SIV-0085 | Syringe Luer Lok 50 mL | | | 3 | | | |

OSTOMY SUPPLIES & ACCESSORIES

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|-------------|------------|----------------------|------|-----|-------------------|-------|
|---------------------|-------------|------------|----------------------|------|-----|-------------------|-------|

OSTOMY SUPPLIES & ACCESSORIES: Products such as ostomy pouches, skin barriers and accessories used in the care & management of urostomies, ileostomies and colostomies. Authorization for Home and Community Care Support Services North East-funded supplies as per approved Ostomy policy and only for short-term treatment to allow patient some time to coordinate private pay options, insurance coverage and ADP applications. Coloplast is the primary vendor. Ostomy supplies required for long-term fistula wound care managment may exceed maximums.

Ostomy Skin Barriers for 2 Piece System: Cut-to-fit flanges (also called "basenlates"), that attach directly to the skin. Measure the stoma accurately and select a cut-to-fit harrier opening which fits

| Connoctably | around the stoma without actually touching it. Too | ` | | | 1 | 1 | | |
|-------------|--|-------------|-----------|-----------------|-------|-------------|--|-----|
| SOS-0068 | 2 Piece Click Barrier Green Flat Cut-to-Fit | SenSura Mio | Coloplast | Coupling/Flange | Υ | | Non-Convexity, Stoma size 10-35mm (3/8"- 1 3/8") | |
| | 21 1000 Click Barrier Green Flat Cut to Fit | Click | 10502 | 40 mm (1 9/16") | | | Tron convexity, eterna size to commit (e/e 1 e/e) | |
| SOS-0025 | 2 Piece Click Barrier Red Flat Cut-to-Fit | SenSura Mio | Coloplast | Coupling/Flange | Υ | | Non-Convexity, Stoma size 10-45mm (3/8"- 1 3/4") | |
| | | Click | 10512 | 50 mm (2") | | | Tron comonity, eterna ella re remin (e/e r e/r / | 9 |
| SOS-0026 | 2 Piece Click Barrier Yellow Flat Cut-to-Fit | SenSura Mio | Coloplast | Coupling/Flange | Υ | | Non-Convexity, Stoma size 10-65mm (3/8"- 2 1/5") | |
| | | Click | 10531 | 70 mm (2 3/4") | | _ | The company, etc. a committee of the com | |
| SOS-0027 | 2 Piece Click Barrier Green Convex Light Cut-to- | SenSura Mio | Coloplast | Coupling/Flange | Υ | | Light Convexity, Stoma size 10-20mm (3/8"- 3/4") | |
| | Fit | Click | 16901 | 40 mm (1 9/16") | | | Light Convexity, Clotha Cl20 To 2011111 (c/c - c/ 1) | |
| SOS-0028 | 2 Piece Click Barrier Red Convex Light Cut-to-Fit | SenSura Mio | Coloplast | Coupling/Flange | V | | Light Convexity, Stoma size 15-30mm (5/8"- 1 1/8") | |
| 000 0020 | · | Click | 16911 | 50 mm (2") | | | Eight Convexity, Cloma Size 10 Comm (6/6 1 1/6) | |
| SOS-0029 | 2 Piece Click Barrier Blue Convex Light Cut-to- | SenSura Mio | Coloplast | Coupling/Flange | V | | Light Convexity, Stoma size 15-40mm (5/8"- 1 9/16") | |
| 000 0023 | Fit | Click | 16921 | 60 mm (2 3/8") | ' | 5 / Week | Light Convexity, Stoma Size 13 40mm (3/6 1 3/10) | |
| SOS-0030 | 2 Piece Click Barrier Yellow Convex Light Cut-to- | | Coloplast | Coupling/Flange | V | | Light Convexity, Stoma size 15-53mm (5/8"- 2 1/16") | |
| 000 0000 | Fit | Click | 16931 | 70 mm (2 3/4") | ' | Maximum 10 | Light Convexity, Stoma Size 13 33mm (3/6 2 1/10) | |
| SOS-0031 | 2 Piece Click Barrier Green Convex Deep Cut-to- | SenSura Mio | Coloplast | Coupling/Flange | | per patient | Deep Convexity, Stoma size 10-20mm (3/8"- 3/4") | |
| 303-0031 | Fit | Click | 16941 | 40 mm (1 9/16") | I | admission | Deep Convexity, Stoma Size 10-20mm (5/6 - 5/4) | |
| SOS-0032 | 2 Piece Click Barrier Red Convex Deep Cut-to- | SenSura Mio | Coloplast | Coupling/Flange | | of any | Deep Convexity, Stoma size 15-30mm (5/8"- 1 1/8") | |
| 303-0032 | Fit | Click | 16951 | 50 mm (2") | ı | combination | Deep Convexity, Stoma Size 13-30mm (3/8 - 1 1/8) | |
| SOS-0033 | 2 Piece Click Barrier BlueConvex Deep Cut-to- | SenSura Mio | Coloplast | Coupling/Flange | ٧ | of sizes | Deep Convexity, Stoma size 15-40mm (5/8"- 1 9/16") | |
| 303-0033 | Fit | Click | 16961 | 60 mm (2 3/8") | Ť | | Deep Convexity, Stoma Size 15-40mm (5/6 - 1 9/16) | |
| SOS-0034 | 2 Piece Click Barrier Yellow Convex Deep Cut-to- | SenSura Mio | Coloplast | Coupling/Flange | γ | | Doon Conveyity Stome size 15 52mm (5/9" 2 1/16") | |
| 303-0034 | Fit | Click | 16971 | 70 mm (2 3/4") | ľ | | Deep Convexity, Stoma size 15-53mm (5/8"- 2 1/16") | |
| 000 0017 | O Diseas Hallisten Dennien Coreas Flat Cort to fit | Carra Divis | Hollister | Coupling/Flange | V | | Non-Convexity, Stoma size up to 32mm (1 1/4") | |
| SOS-0017 | 2 Piece Hollister Barrier Green Flat Cut-to-fit | CeraPlus | 11202 | 44 mm (1 3/4") | ľ | | Infused with ceramide to protect peristomal skin | |
| 000 0010 | O Pierra Hallistan Bandan Bad Flat O at ta St | O Dl | Hollister | Coupling/Flange | V | | Non-Convexity, Stoma size up to 44mm (1 3/4") | 0/2 |
| SOS-0018 | 2 Piece Hollister Barrier Red Flat Cut-to-fit | CeraPlus | 11203 | 57 mm (2 1/4") | Y | | Infused with ceramide to protect peristomal skin | |
| 000 0010 | O. Pierra Hallistan Bandan Orana Orana Orana St. | O DI | Hollister | Coupling/Flange | \ \ \ | 1 | Convex, Stoma size up to 25mm (1") | |
| SOS-0019 | 2 Piece Hollister Barrier Green Convex Cut-to-fit | CeraPius | 11402 | 44 mm (1 3/4") | Υ | | Infused with ceramide to protect peristomal skin | |
| 000 0000 | O Pierra Hallistan Bandan Bad Oannan C. 11. (i) | O DI | Hollister | Coupling/Flange | V | 1 | Convex, Stoma size up to 38mm (1 1/2") | |
| SOS-0020 | 2 Piece Hollister Barrier Red Convex Cut-to-fit | CeraPlus | 11403 | 57 mm (2 1/4") | Y | | Infused with ceramide to protect peristomal skin | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes | | | | | |
|---------------------|--|---------------------------|----------------------|-----------------------------------|-----|-----------------------------|--|--|--|--|--|--|
| | rostomy Pouches for 2 Piece System: Order with appropriate brand and size Ostomy Skin Barrier A drainable pouch that connects to the barrier flange to collect urine with a integrated drain valve at the bottom to empty the bag. | | | | | | | | | | | |
| SOS-0035 | 2 Piece SenSura Mio Click Urostomy Pouch Green | SenSura Mio Click Maxi | Coloplast 11492 | Coupling/Flange 40mm (1 9/16") | Υ | 5 / Week | Bag length 26cm (10 1/4"), | | | | | |
| SOS-0036 | 2 Piece SenSura Mio Click Urostomy Pouch Red | SenSura Mio Click Maxi | Coloplast 11498 | Coupling/Flange 50mm (2") | Υ | Maximum 10 per patient | Transparent | | | | | |
| SOS-0021 | 2 Piece Hollister Urostomy Pouch Green | Hollister | Hollister 18902 | Coupling/Flange 44mm (1 3/4") | Υ | admission of any | | | | | | |
| SOS-0022 | 2 Piece Hollister Urostomy Pouch Red | Hollister | Hollister 18903 | Coupling/Flange 57mm (2 1/4") | Υ | combination of sizes | | | | | | |
| | Pouches for 2 Piece System: Order with a sto the barrier flange. An integrated closure or class | | | | | ouch with a larg | ger opening for ileostomies and colostomies | | | | | |
| SOS-0037 | 2 Piece SenSura Mio Click Drainable Pouch Green | SenSura Mio Click Maxi | Coloplast 11463 | Coupling/Flange 40mm (1 9/16") | Υ | | | | | | | |
| SOS-0038 | 2 Piece SenSura Mio Click Drainable Pouch Red | SenSura Mio Click Maxi | Coloplast 11473 | Coupling/Flange 50mm (2") | Υ | 5 / Week | Bag length 28cm (11") Transparent | | | | | |
| SOS-0039 | 2 Piece SenSura Mio Click Drainable Pouch Blue | SenSura Mio Click Maxi | Coloplast 11483 | Coupling/Flange 60mm (2 3/8") | Υ | Maximum 10 | Easiclose™ Wide Outlet | | | | | |
| SOS-0040 | 2 Piece SenSura Mio Click Drainable Pouch Yellow | SenSura Mio Click Maxi | Coloplast 11487 | Coupling/Flange 70mm (2 3/4") | Υ | per patient admission | | | | | | |
| SOS-0023 | 2 Piece Hollister Drainable Pouch Green | New Image Lock N Roll | Hollister 18132 | Coupling/Flange 44mm (1 3/4") | Υ | of any combination of sizes | May only be ordered with the Hollister 2-Piece Barrier Bag length 30cm (12") | | | | | |
| SOS-0024 | 2 Piece Hollister Drainable Pouch Red | New Image Lock N Roll | Hollister 18133 | Coupling/Flange 57mm (2 1/4") | Y | | Transparent No filter | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes | | | | | |
|----------------------------|--|-----------------|----------------------|--|----------------|------------------------------------|---|------|--|--|--|--|
| One Piece | Piece Urostomy System: Includes the cut-to-fit ostomy skin barrier and urostomy pouch as a single unit to collect urine with an integrated drain valve at the bottom to empty the bag. | | | | | | | | | | | |
| SOS-0041 | 1 Piece Urostomy Pouch Convex Light 10-23mm Cut-to-Fit | SenSura Mio | Coloplast 16835 | Stoma Size 10-23mm (3/8-7/8") | Υ | | 1 | | | | | |
| SOS-0042 | 1 Piece Urostomy Pouch Convex Light 15-33mm Cut-to-Fit | SenSura Mio | Coloplast 16836 | Stoma Size 15-33mm (5/8-1 1/4") | Υ | 5 / Week | Light Convexity Bag length 26cm (10 1/4") Transparent | | | | | |
| SOS-0043 | 1 Piece Urostomy Pouch Convex Light 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16837 | Stoma Size 10-43mm (3/8-1 3/4") | Υ | Maximum 10 per patient | Transparent | | | | | |
| SOS-0044 | 1 Piece Urostomy Pouch Convex Deep 10-23mm Cut-to-Fit | | Coloplast 16865 | Stoma Size 10-23mm (3/8-7/8") | Υ | admission of any | David Committee | 0.0) | | | | |
| SOS-0045 | 1 Piece Urostomy Pouch Convex Deep 15-33mm Cut-to-Fit | | Coloplast 16866 | Stoma Size 15-33mm (5/8-1 1/4") | Υ | combination of sizes | Deep Convexity Bag length 26cm (10 1/4") Transparent | YS7 | | | | |
| SOS-0046 | 1 Piece Urostomy Pouch Convex Deep 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16867 | Stoma Size 10-43mm (3/8-1 3/4") | Υ | | Transparent | | | | | |
| One Piece until it require | • | able pouch with | a larger opening fo | or ileostomies and colostor | mies. <i>A</i> | n integrated ro | oll closure or clamp is used to keep the pouch closed | | | | | |
| SOS-0047 | 1 Piece Drainable Pouch Flat 10-55mm Cut-to- Fit | SenSura Mio | Coloplast 10481 | Stoma Size 10-55mm (3/8-1 1/8") | Υ | | | | | | | |
| SOS-0048 | 1 Piece Drainable Pouch Convex Light 10-23mm Cut-to-Fit | SenSura Mio | Coloplast 16737 | Stoma Size 10-23mm (3/8-7/8") | Υ | | | | | | | |
| SOS-0049 | 1 Piece Drainable Pouch Convex Light 15-33mm Cut-to-Fit | SenSura Mio | Coloplast 16738 | Stoma Size 15-33mm (5/8-1 1/4") | Υ | 5 / Week | Bag length 28cm (11") | | | | | |
| SOS-0050 | Cut-to-Fit | SenSura Mio | Coloplast 16739 | Stoma Size 10-43mm (3/8-1 3/4") | Υ | Maximum 10 | | | | | | |
| SOS-0051 | 1 Piece Drainable Pouch Convex Deep 10-23mm Cut-to-Fit | | Coloplast 16755 | Stoma Size 10-23mm (3/8-7/8") | Υ | per patient admission of any | Transparent or Opaque (with inspection window) | | | | | |
| SOS-0052 | 1 Piece Drainable Pouch Convex Deep 15-33mm Cut-to-Fit | SenSura Mio | Coloplast 16756 | Stoma Size 15-33mm (5/8-1 1/4") | Υ | combination of sizes | | | | | | |
| SOS-0053 | 1 Piece Drainable Pouch Convex Deep 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16767 | Stoma Size 10-43mm (3/8- 1 3/4") | Υ | | | | | | | |
| SOS-0071 | 1 Piece High Output Pouch 10-100mm Cut-to-Fit | SenSura Mio | Coloplast 18666 | Baseplate hole size: 10-100mm - 725ml | Υ | | Trasparent cut-to-fit 10-100mm | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | | | | | | |
|---------------------|---|-------------|----------------------|--------------------------------|-----|----------------------------------|--|--|--|--|--|--|--|
| | ediatric 2 Piece System: Two-piece system with adhesive coupling. Soft and comfortable - as soft and flexible as a one-piece pouch, ensuring freedom of movement while still wearing a two-piece. For childrenged aged six months to four years. The pediatric system is simple and easy to use and the coupling flange can be lifted so that the pouch can be applied without pressure on the peristomal area. | | | | | | | | | | | | |
| SOS-0069 | Kids Flex Barrier Flat Cut-to-Fit | SenSura Mio | Colplast 18747 | Baseplate hole size: 0-35mm | Υ | 5 / Week Maximum 10 per patient | Cut- to- fit 0-35 mm Designed for kids aged six months to four years, to be worn in combination with a compatible SenSura Mio Kids Drainable Pouch SSO-0070 | | | | | | |
| SOS-0070 | Kids Flex Drainable Pouch Midi | SenSura Mio | Coloplast 18724 | Coupling Size 35 mm | Y | admission of any | Opaque 35 mm Designed for kids aged six months to four years, to be worn in combination with a compatible SenSura Mio Kids barrier SSO-0069 | | | | | | |
| Ostomy Ad | ccessories: | | | | | | | | | | | | |
| SOS-0054 | Ostomy Support Belt for SenSura Mio Standard | Brava Belt | Coloplast 0423 | Each up to 100 cm (40") | Υ | 1 per Patient | Standard Size | | | | | | |
| SOS-0055 | Ostomy Support Belt for SenSura Mio XXL | Brava Belt | Coloplast 0424 | Each up to 135 cm (61") | Υ | Admission | XXL Size | 7 | | | | | |
| SOS-0056 | Ostomy Strip Paste | Brava | Coloplast 2655 | Each 57g (2 oz) | Υ | 5 / Patient Admission | Moldable. Fills gaps to create a tight seal to prevent leakage. | Brave* Strip Partir & Complete | | | | | |
| SOS-0057 | Ostomy Protective Sheet | Brava | Coloplast 3210 | Each 10 x10 cm 4 x 4" | Υ | 3 | Absorbs moisture and protects peristomal skin when used under the barrier to minimize risk of maceration | Breve* Productive Short © Coloptast | | | | | |
| SOS-0058 | Ostomy Paste (Tube) | Brava | Coloplast 12050 | Each Tube 60g (2 oz) | Υ | 1 / Patient Admission | Alcohol-free paste fills folds and creases to maintain a tight fit between the stoma and the ostomy barrier to reduce leakage | Minore Particular Part | | | | | |
| SOS-0059 | Ostomy Powder | Brava | Coloplast 1907 | Each Bottle 25g (1 oz) | Υ | 1 / Patient Admission | It absorbs moisture so the skin is kept dry to reduce irritation. Can be used around the stoma or under the barrier baseplate. | Report Poorder Poorder Octobross | | | | | |
| SOS-0066 | Mouldable Skin Barrier Ring | Brava | Coloplast 12045 | 18 mm diameter 4.2 mm thick | Y | 1 | A durable ring for reduced leakage. Creates a durable seal between the stoma and the base plate, thus reducing leakage. It is easy to use and can be molded to fit the patient's body shape. | 0 | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|-------------------|-----------------------|---|-----|-------------------|--|
| SOS-0065 | Ostomy Care Kit | Coloplast | Coloplast | Each Kit | Υ | 1 | Maximum Quantity: 1 Kit/Patient Admission |
| One piece | e ostomy flat system | | | 1 | | | E RIEMEN A |
| One piece | e ostomy convx light system | | | 1 | | | Colopiast |
| Scissors | | | | 1 | | | |
| Mirror | | | | 1 | | | NAME OF THE PARTY |
| Adhesive | Remover Wipe | | | 1 | | | Allerent Control of the Control of t |
| Lubricatin | g Deordorant | | | 1 | | | |
| Measurin | g Guide & Marker | | | 1 | | | Stomaguide |
| Coloplast | Care Information and Forms (Billingual) | | | 1 | | | Coroplast |
| SOS-0060 | Fistula and Wound Pouch | Eakin | ConvaTec 839264 | For wounds up to 245 x 160 mm (9.7" x 6.3") | Υ | 5 | Flexible pouches designed for skin protection and contain drainage from wounds, fistulas and difficult stoma sites. Horizontal cutting area. Suitable with remote drainage attachment, tap closure & window for wound access. |
| SOS-0061 | Skin Barrier Ring Small | Eakin Cohesive | ConvaTec 839002 | Each 50 mm (2") | Υ | 1 | Can be moulded to a variety of shapes to help prevent leakage around the stoma by acting as filler |
| SOS-0062 | Skin Barrier Ring Large | Eakin Cohesive | ConvaTec 839001 | Each 100 mm (4") | Υ | 1 | for uneven skin surfaces. |
| SOS-0063 | Antimicrobial Ostomy Ring | Hydrofera Blue | Hydrofera HBRF2650 | 60 mm (2.5") Diameter | Υ | 3 | - Not to be used longer than 6 months Not to be use on pregnant or nursing patients. |
| SOS-0064 | Ceramide Skin Barrier Ring | Adapt CeraRing | Hollister 8805 | 48 mm (2") Diameter | Υ | 1 | Infused with ceramide to help protect the peristomal skin's natural moisture barrier to prevent leakage and dryness. |

URINARY & DRAINAGE SUPPLIES

| CHRIS | Description | Brand Name | Manufacturer | Size | APO | Weekly | Notes |
|------------|-------------|------------|--------------|------|-----|---------|-------|
| Order Code | _ | | Code | | | Maximum | |

Urinary Supplies: Medical supplies that helps drain urine from the bladder. It is used by some patients who have had urological or gynecological surgery, or who have a condition that makes urination difficult. Home and Community Care Support Servcies North East-funded supplies may be authorized only while the patient is in need of nursing services. ODSP, social assistance, private insurance, DVA must be accessed as possible sources of funding on a long term basis as per approved policy. Most commonly used sizes are listed below. All other sizes (including pediatric) and coatings are available for special order only as per the Exception Item process.

FOLEY CATHETERS: 2-way Foley (inflation and drainage) indwelling catheters. For short- or long-term catheter patients, the catheter should be the smallest size diameter that permits good urine flow. If a larger diameter was used to facilitate passage of blood clots, the size can be decreased by 2-4 Fr per catheter change until the optimal size is reached. Likewise, balloon sizes should be small, with 5 mL balloons, reserving the 30 mL only for bleeding episodes, usually postoperatively. Long term catheters (over 1 months use) often changed every 4 weeks. Patients with frequent blockage may need catheter changed every 2-3 weeks or more often.

| SUR-0001 | Foley 2 Way Hydrogel Coated Latex 12Fr 5 mL | BARD Canada | BARD 226512 or equivalent | 12 Fr 5 mL Balloon | Υ | 2/month | |
|----------|---|-------------|---------------------------------|------------------------|---|---------|--|
| SUR-0002 | Foley 2 Way Hydrogel Coated Latex 14Fr 5 mL | BARD Canada | BARD 226514 or equivalent | 14 Fr 5 mL Balloon | Υ | 2/month | |
| SUR-0003 | Foley 2 Way Hydrogel Coated Latex 16Fr 5mL | BARD Canada | BARD 226516 or equivalent | 16 Fr 5 mL Balloon | Υ | 2/month | |
| SUR-0004 | Foley 2 Way Hydrogel Coated Latex 18Fr 5 mL | BARD Canada | BARD 226518 or equivalent | 18 Fr 5 mL Balloon | Υ | 2/month | |
| SUR-0014 | Foley 2 Way Hydrogel Coated Latex 14Fr 30mL | BARD Canada | BARD 226614 or equivalent | 14 Fr 30 mL Balloon | Υ | 2/month | |
| SUR-0015 | Foley 2 Way Hydrogel Coated Latex 16Fr 30 mL | BARD Canada | BARD 226616 or equivalent | 16 Fr 30 mL Balloon | Υ | 2/month | |
| SUR-0016 | Foley 2 Way Hydrogel Coated Latex 18Fr 30 mL | BARD Canada | BARD 226618 or or equivalent | 18 Fr 30 mL Balloon | Υ | 2/month | |
| SUR-0017 | Foley 2 Way Hydrogel Coated Latex 20Fr 30 mL | BARD Canada | BARD 226620 or equivalent | 20 Fr 30 mL Balloon | Υ | 2/month | |
| SUR-0116 | Foley 2 Way Red Latex Tiemann Coude Tip 14Fr 5 mL | BARD Canada | BARD 0102L14 | 14 Fr 5 mL Balloon | Υ | 2/month | For difficult insertions, coude-tip Foley catheters feature a bent tip that maneuvers around obstructions leading to the bladder (eg. men with |
| SUR-0117 | Foley 2 Way Red Latex Tiemann Coude Tip 16Fr 5 mL | BARD Canada | BARD 0102L16 | 16 Fr 5 mL Balloon | Υ | 2/month | enlarged prostate) |
| SUR-0027 | Foley 2 Way 100% Silicone 12 Fr 10mL | Folysil | Coloplast AA6112 | 12 Fr 10 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0028 | Foley 2 Way 100% Silicone 14Fr 10mL | Folysil | Coloplast AA6114 | 14 Fr 10 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0037 | Foley 2 Way 100% Silicone 16Fr 15mL | Folysil | Coloplast AA6116 | 16 Fr 15 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0046 | Foley 2 Way 100% Silicone 18Fr 15mL | Folysil | Coloplast AA6118 | 18 Fr 15 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0033 | Foley 2 Way 100% Silicone 18Fr 30mL | Folysil | Coloplast AA6C18 | 18 Fr 30 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|---|---------------------------|------------------------|----------------------------|-----------|-------------------|--|
| SUR-0034 | Foley 2 Way 100% Silicone 20Fr 30mL | Folysil | Coloplast AA6C20 | 20 Fr 30 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0047 | Foley 2 Way 100% Silicone Coude Tip 14Fr 15mL | Folysil | Coloplast AA6314 | 14 Fr 15 mL Balloon | Υ | 2/month | Tiemann (coude) tip for difficult insertions. Features a bent tip that maneuvers around obstructions leading to the bladder (eg. men with |
| SUR-0121 | Foley 2 Way 100% Silicone Coude Tip 16Fr 15mL | Folysil | Coloplast AA6316 | 16 Fr 15 mL Balloon | Υ | 2/month | enlarged prostate). Latex free |
| | ENT CATHETERS: Closed system kit designed for stern tract infections. | erile intermittent o | atheterization. All-in | n-one ready to use se | elf-lubic | ating cathete | er and sterile bag set. The unique closed system helps to reduce the |
| SUR-0112 | Female Intermittent Catheter 10 Fr | SpeediCath Compact Set | Coloplast 28520 | 10 Fr x 9 cm | Υ | 35 | |
| SUR-0113 | Female Intermittent Catheter 12 Fr | SpeediCath Compact Set | Coloplast 28522 | 12 Fr x 9 cm | Υ | | Bag is directly attached to the catheter as a closed system Single patient use - Disposable. |
| SUR-0114 | Female Intermittent Catheter 14 Fr | SpeediCath Compact Set | Coloplast 28524 | 14 Fr x 9 cm | Υ | 35 | Available in male and female versions Maxiumum 5/day |
| SUR-0115 | Male Intermittent Catheter 12-18 Fr | SpeediCath Compact Set | Coloplast 28422 | 12-18 Fr 30 cm long | Υ | 35 | |
| SUR-0125 | Male External Condom Catheter Silicone 23mm | Freedom Clear | Coloplast 505711 | 23 mm | Υ | 7 | Standard Adhesive; Non-sterile; Single Patient Use; Disposable |
| SUR-0126 | Male External Condom Catheter Silicone 28mm | Freedom Clear | Coloplast 505731 | 28 mm | Υ | 7 | Standard Adhesive; Non-sterile; Single Patient Use; Disposable |
| SUR-0127 | Male External Condom Catheter Silicone 31mm | Freedom Clear | Coloplast 505741 | 31 mm | Υ | 7 | Standard Adhesive; Non-sterile; Single Patient Use; Disposable |
| | PRAINAGE BAGS: Urine drainage bags are used to wist valve) are available by an approved Exception | | • | • . | | | for use with Foley Catheters. Drainage bags with slide valve |
| SUR-0085 | Leg Bag Kit 600 mL Flip-Flo Valve with Straps and 18" Tubing | Bard Canada | Dispoz-A-Bag 150180 | 19oz / 562ml 18" Tubing | Υ | 2 | - For long-term catheter users, the bags should be replaced every 5-7 days |
| SUR-0086 | Leg Bag 562 mL Flip-Flo Valve with Straps | Bard Canada | Dispoz-A-Bag 1501MF | 19 oz / 562ml Each | Υ | 2 | - Bags must be emptied when 1/3 to 1/2 full - Adjustable straps are included and attached to the leg bags. |
| SUR-0089 | Leg Bag 946mL Flip Flo Valve with Straps | Bard Canada | Dispoz-A-Bag 1501LF | 32 oz / 946ml Each | Υ | 2 | No longer best practice to change bag for overnight use unless the volume requires this to prevent leakage/overflow; |
| SUR-0090 | Drainage Bag 2000mL with Anti-Reflux Device, Catheter tubing and Metal Drain Clamp | Bard Canada | Bard 154002 | 2000 mL Each | Υ | 2 | - Larger sized 2000 ml bag often used for patients who are bedridden or end of life. |
| SUR-0094 | Urinary Leg Bag Extension Tubing 18" | Medline | Medline DYND12550 | 18 inches Each | Υ | 2/month | 18" extension tubing can be cut to preferred length to meet patient needs |
| SOS-0013 | Nephrostomy Drainage Bag 600 mL | Argon | DBAG600 | 600 mL | Υ | 1 | With adsjustable elastic strap for waist or leg. Provide 2 at admission then 1 / week until discharge. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|-------------------------|--------------------------------|---|-------|-------------------|--|
| CATHETER | STABILIZATION: Adhesive devices used to stabilize | ze indwelling ur | inary catheters, re | ducing catheter mov | vemer | nt and minin | nizing accidental dislodgments |
| SUR-0111 | Universal Foley Catheter Stabilization Device | StatLock | Bard FOL0102 | Each | Υ | 2/month | Maximum for nephrostomy care is up to 4/mth |
| SUR-0104 | Statlock Universal Plus Catheter Stabilization Device Large 12-14 FR | StatLock | Bard VUPD1214 | 12 - 14 fr | Υ | 2/month | Destricts and between the grad file viet valle all and between airca |
| SUR-0105 | StatLock Universal Plus Catheter Stabilization Device Medium 10 - 12 FR | StatLock | Bard VUPD1012 | 10 - 12 fr | Υ | 2/month | Restricts catheter motion and fits virtually all catheter sizes. For biliary and nephrostomy drainage and peritoneal lavage catheters. |
| SUR-0106 | StatLock Universal Plus Catheter Stabilization Device Small 6 - 8.5 FR | StatLock | Bard VUPD68 | 6 - 8.5 fr | Υ | 2/month | |
| CATHETER | ACCESSORIES: There are various accessories the | at can be used | to assist with eithe | r indwelling cathete | r use | or the proce | ess of intermittent catheterisation. |
| SUR-0078 | Catheter Plug and Drain Tube Cover Set | Medline | Medline DYND12200 | Each | Υ | 2/month | Fits securely into funnel end of most Foley catheters; Sterile; Cover protects connector on drainage bag |
| SUR-0110 | Syringe Catheter Tip 50 mL | BD | BD 309620 | 50 mL | Υ | 7 | Used for flushing/cleaning catheters, gastronomy tubes and other devices. |
| SUR-0120 | Bulb Irrigation Syringe 60 mL | Medline | Medline DYND20125 | 60 mL | Υ | 7 | Soft, pliable bulb for maximum control of drawn or expelled volume. Easy to grip and use with only one hand; Sterile |
| SUR-0091 | Sterile Bacteriostatic Lubricant Water Based | Healthcare Plus | Healthcare Plus 248-LUB-035 | 3.5 g pouch | Υ | 2/month | For use only with Foley catheters. Not required for SpeediCath lubricated intermittent catherters. |
| SUR-0103 | Drainage Tube Universal 1 Way Stopcock | MED-Rx | 10-1210 | 16 inch | Υ | 2/month | Maximum for nephrostomy care is up to 4/mth |
| 3011-0103 | Brainage Tube Oniversal T way Stopcock | Cook Medical | G02278 | 30 cm | Υ | 2/month | Due to recurring backorders this is a substitute for MED-RX 10-1210 |
| | Specimen container | Star Plex Scientific | Fisher Scientific 360993197 | 90 mL | Υ | 2 | Sterile polypropylene container for collection, transport and storage of specimens, Single use, Disposable |
| | US BLADDER IRRIGATION (CBI): Approved for tin nent may be as frequent at 1 bag/hour per day. | ne-limited treati | ment for Continuou | ıs Bladder Irrigation | only. | Usually red | uired to discharge patient from hospital for end of life or hospice |
| SUR-0122 | 0.9% Sodium Chloride Irrigation Solution 3000 mL | Baxter | Baxter JB7127I | 3000 mL | N | by | Must be authorized/ordered by Care Coordinator SPECIAL ORDER ITEM: Availability of stock with the vendor needs to be confirmed prior to ordering. |
| SUR-0123 | Bladder Irrigation Tubing Set | Baxter | Baxter JC4040 | 0.195" Inner Dia. Length: 2.1m (81" Long) | Υ | 1 | Includes tubing, drip chamber, regulating clamp, and catheter adapter. Plug directly into 3-way catheter or us catheter adaptor unit when 2-way catheter is used. |
| SUR-0124 | Foley 3 Way Catheter for CBI Silicone 22 Fr | X-Flow | Coloplast AB6122 | 22 Fr 30-50 mL Balloon | Υ | 1 | Patient usually discharged from hospital with required catheter. Couvalier tip appropriate for both male and female patients. Additional patient-specific sizes available with Exception Item Request approval |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|---|-----------------|----------------------|----------------------|---------|-------------------|--|
| | ITS: Disposable, Single use prepared kits of comm Additional items may also be ordered as needed. | only used items | | ich kit can be order | ed indi | | n entire kit is not required if an individual CHRIS order code is |
| SUR-0079 | Foley Catheter Tray | MED-Rx | CHS 85-5050 | Each Tray | Υ | 2 / Month | For catherter insertion 1/month |
| | CSR Wrap Blue |] | | 1 | | | Order catheters, leg bags or drainage bags seperately |
| | Plastic Transfer Forcep | | | 1 | | | (may require 1 - SIV-0082 - 10mL Luer Lok syringe to deflate balloon for removal of catheter) |
| | Graduated Basin | | | 1 | | | If irrigation required add: |
| | Gloves Medium (Pair) Non-latex | | | 1 | | | SUR-0080 Irrigation Tray Kit with 60 mL Piston Syringe and SSO-0019 Normal Saline 500mL bottle |
| | Fenestrated Drape Paper | | | 1 | | | |
| | Towel Paper | | | 1 | | | |
| | Povidone Iodine Swabstick (3 pack) |] | | 2 | | | |
| | Lubricant Water Based |] | | 1 | | | |
| | Foley Catheter Stabilization Device |] | | 1 | | | |
| | Sterile Water 10 mL Pre-Filled Syringe |] | | 1 | | | |
| SUR-0080 | Irrigation Tray Kit with 60 mL Piston Syringe | MED-Rx | CHS 85-5042 | Each Kit | Υ | 1 | Disposable, single use |
| | Irrigation Tray | | | 1 | | | |
| | 500 mL Graduated Container |] | | 1 | | | |
| | 60 mL Piston Irrigation Syringe with Cap |] | | 1 | | | |
| | Alcohol Prep Pad |] | | 1 | | | |
| | PVP Prep Pad |] | | 1 | | | |
| | Underpad 17" x 18" | | | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|------------|-------------------------|----------|-----|-------------------|--|
| SMK-0008 Pe | MK-0008 Percutaneous Nephrostomy Kit: | | Various - See each item | Each Kit | Υ | 1 | Patient usually discharged from hospital with required catheter. Additional patient-specific sizes available as per this catalog or by |
| SDR-0005 | Dressing Tray Kit - Disposable Single Use | | | 2 | | | Exception Item Request approval |
| SDR-0055 | Drain Gauze Sponge Trach Sterile Woven 4"x4" | | | 2 | | | |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | | | 2 | | | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | | | 2 | | | |
| SSO-0021 | Normal Saline 10 mL Pre-Filled Syringe | | | 1 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 14 | | | |
| SUR-0085 | Leg Bag Kit 600 mL Flip-Flo Valve with Straps and 18" Tubing | | | 1 | | | |
| SWC-0128 | Transparent Film Dressing 15x20cm | | | 2 | | | |
| SUR-0103 | Drainage Tube Universal 1 Way Stopcock | | | 2 | | | |
| SMK-0044 Pa | aracentesis Drainage Kit | | | | Υ | 1 | To be ordered by and provided to Nurse Practitioners or Physicians only for abdomen or chest drainage procedures for eligible palliative |
| Not Applicable | Paracentesis Drainage Tray w/9cm Catheter | Smartap | CardioMed CM-0020T | 1 | | | patients to support end of life care. Home and Community Care Support Services North East will provide |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | the medical supplies to the physician for eligible patients only. Funding or supplies for topical or local anesthetic for drainage procedures is not available and is to be provided by the physician or is |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | | | 2 | | | at the cost of the patient. |
| SCM-0019 | Blue Absorbent Pad 23"x24" | | | 2 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 2 | | | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | | | 2 | | | |

FEEDING, SUCTION & TRACHEOSTOMY SUPPLIES

| CHRIS | Description | Brand Name | Manufacturer | Ci=o | 400 | Weekly | Notes |
|------------|-------------|------------|--------------|------|-----|---------|-------|
| Order Code | Description | Brand Name | Code | Size | APO | Maximum | Notes |

Initial supply authorization for four weeks only for short-term treatment and to allow patients to coordinate insurance or ADP applications. Supplies will only continue to be provided as per approved policy. Funded feeding supplies shall only be provided for Kangaroo ePump or Joey Pump and may not be compatible with other pumps. Home and Community Care Support Services North East is not obligated to provide supplies beyond defined periods while patients wait for third-party funding.

ENTERAL FEEDING

| ETM-1002 | | | NOTE: Equipment is available for rental as per approved policies from Medigas (Praixair). Please review to the Medigas Medical Equipment Catalogs for more detailed information before ordering. | | | | | | | | | |
|----------|--|-----------------------|--|--------------------------------|---------|------------------------|---|--|--|--|--|--|
| ETM-1005 | IV pole - portable | Please review t | o the Medigas Med | lical Equipment Cata | alogs f | or more detail | ed information before ordering. | | | | | |
| SFS-0005 | Enteral Feeding Gravity Bag | Kangaroo | Cardinal Health Z8884702500 | 1000 ml | Υ | 7 | To be used in the absence of a feeding pump. Non-sterile | | | | | |
| SFS-0006 | Enteral Feeding Pump Set 500 ml | Kangaroo Joey | Cardinal Health Z762055 | 500 ml | Υ | 7 | | | | | | |
| SFS-0031 | Enteral Feeding Pump Set 1000 ml | Kangaroo Joey | Cardinal Health Z763656 | 1000 ml | Υ | 7 | | | | | | |
| SFS-0035 | Enteral Feeding Pump Set with Flush Set 1000ml | Kangaroo Joey | Cardinal Health Z763662 | 1000 ml | Υ | 7 | For continuous overnight feeds ONLY | | | | | |
| SFS-0045 | Feeding Tubing Set with Safety Screw Spike | Kangaroo Joey | Cardinal Health Z765559 | 1 Tube | Υ | 7 | By request from Physician or Registed Dietician only for patients with a clinical need for a closed system. Tubing only - no bag. | | | | | |
| SFS-0010 | Gastrostomy Tube with Y Port 16 Fr | Kangaroo | Cardinal Health Z8884720163E | 16 Fr x 9" 20 cc | Υ | 1/Patient Admission | | | | | | |
| SFS-0036 | Gastrostomy Tube with Y Port 20 Fr | Kangaroo | Cardinal Health Z8884720205E | 20 Fr x 9" 20 cc | Υ | 1/Patient Admission | | | | | | |
| SFS-0012 | Nasogastric Feeding Tube 5 Fr (Infant) | Kangaroo | Cardinal Health Z461412E | 5 Fr x 20" 1.7 mm x 51 cm | Υ | 1/month | Remains in situ for 30 days Purple Polyurethane Sterile | | | | | |
| SFS-0037 | Nasogastric Feeding Tube PVC 5 Fr (Infant) | Kangaroo | Cardinal Health Z460802E | 5 Fr x 16" 1.7 mm x 41 cm | Υ | 2/week | Used only for patients who pull them out frequently PVC Sterile | | | | | |
| SFS-0038 | Nasogastric Feeding Tube 6.5 Fr (Infant) | Kangaroo | Cardinal Health Z461438E | 6.5 Fr x 36" 2.2 mm x 91 cm | Υ | 1/month | Remains in situ for 30 days Purple Polyurethane Sterile | | | | | |
| SFS-0015 | INISCONSCINC FARMING LUNG 12 Fr | Kangaroo Entriflex | Cardinal Health Z8884721252E | 12 Fr x 43" 4 mm x 109 cm | Υ | 1/month | With stylet and rigid outlet port | | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | | | |
|---------------------|--|--|---------------------------------|--------------------------|-----|------------------------|---|--|--|--|
| SFS-0040 | Adapter Y-Port-G-tube/P.E.G. 16 Fr | Kangaroo | Cardinal Health Z8884751622E | Fits 16 Fr or smaller | Υ | 1/Patient Admission | Non-Sterile. For use with the Kangaroo P.E.G. Kits and can be used as the replacement Y-port on all Kangaroo gastrostomy tubes. 1 per | | | |
| SFS-0041 | Adapter Y-Port-G-tube/P.E.G. 20 Fr | Kangaroo | Cardinal Health Z8884752042E | Fits all 18 Fr or bigger | Υ | 1/Patient Admission | patient admission. Replace only if port is clogged or broken when patient still eligible for supplies. | | | |
| SFS-0044 | Feeding Tube Y-Site Extension Set | Kangaroo | Cardinal Health Z8884705008 | 12 inch | Υ | 1 | Not compatible with low-profile tubes like MIC-KEY or Nutriport | | | |
| SFS-0042 | Purple Enteral Oral Syringe 3 mL | Monoject | Cardinal Health Z403S | 3 mL | Υ | 7 | Standard Tip; Non sterile; Can be reused if washed with same patient For Neonatal and Pediatric Feeding Systems | | | |
| SFS-0043 | Purple Enteral Oral Syringe 10 mL | Monoject | Cardinal Health Z410S | 10 mL | Υ | 7 | Standard Tip; Non sterile; Can be reused if washed with same patient For Neonatal and Pediatric Feeding Systems | | | |
| SFS-0047 | Purple Enteral Oral Syringe 35 mL Sterile with Tip Cap | Monoject | Cardinal Health Z435SG | 35 mL | Υ | 7 | Standard Tip; Sterile; For Neonatal and Pediatric Feeding Systems Single Patient Use with Nasogastric Feeding Tube | | | |
| SFS-0048 | 3 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z403SE | 3 mL | Υ | 7 | Can be reused if washed with same patient ENFit® Connection | | | |
| SFS-0049 | 12 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z412SE | 12 mL | Υ | 7 | Can be reused if washed with same patient ENFit® Connection | | | |
| SFS-0050 | 35 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z435SE | 35 mL | Υ | 7 | Can be reused if washed with same patient ENFit® Connection | | | |
| SFS-0051 | 60 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z460SE | 60 mL | Υ | 7 | Can be reused if washed with same patient ENFit® Connection | | | |
| SFS-0046 | Low Profile Tube Feed Extension Set | MIC-KEY Secur-Lok | Kimberly Clark 0121-12 | 12 inch | Υ | 1 | Primarily used with pediatric patients. SECUR-LOK Right-Angle Connector; 2 Port "Y" connector and Tubing clamp for use with low-profile MIC-KEY tubes. Not compatible with Nutriport low-profile tubes. | | | |
| SFS-0039 | Nasogastric Tube Stabilization Device (Pediatric) | Statlok | Bard NG0100 | Pediatric | Υ | 1 | Order for babies and small children For nasogastric tubes and feeding tubes | | | |
| SFS-0019 | Nasogastric Tube Stabilization Device (Adult) | Statlok | Bard NG0120 | Adult | Υ | 1 | For nasogastric tubes and feeding tubes | | | |
| SUCTION | CATHETERS, TUBING AND ACCESSOR | RIES | | | | | | | | |
| ETM-1006 | Suction machine (includes nebulizer bottle, 2 trach masks, and tubing) | NOTE: Equipment is available for rental as per approved policies from Medigas (Praixair). | | | | | | | | |
| ETM-1015 | Gastric Drainage Aspirator | Please review to the Medigas Medical Equipment Catalogs for more detailed information before ordering. | | | | | | | | |
| SOT-0090 | Disposable Suction Tubing Kit for ETM-1015 | Gomco | Allied 01-90-2000 | 18" and 72" | Υ | 7 | Includes two different length tubes per kit | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes |
|---------------------|--|------------|----------------------|-------------|-----|-------------------|---|
| SOT-0091 | Disposable Filter (Box of 3) for ETM-1015 | Gomco | Allied 01-90-3100 | Box of 3 | Υ | 1 | Filter only needs replacing if it gets wet |
| SOT-0092 | Disposable Collection Bottle for ETM-1015 | Gomco | Allied 01-90-3695 | 1100 mL | Υ | 7 | Plastic; Single Patient Use |
| ISO1-0101 | Connector for Suction Tubing for ETM-1015- ONLY | Medline | DYND50510 | Each | Υ | | 5-in-1 straight plastic connector to connect tubing to ETM-1015 Gastric Drainage Aspirator |
| SOT-0072 | Universal Suction Tubing with Connectors | Ornex | Medline OR66A | 1/4" X 72" | Υ | 1 | Straight connectors feature ribbing which provides a firm connection to suction devices of all kinds. Designed to handle 28 Hg pressure. |
| SOT-0053 | Suction Catheter 8 FR | Med-Rx | CHS 50-1008 | 8 Fr x 15" | Υ | 7 | |
| SOT-0054 | Suction Catheter 12Fr | Med-Rx | CHS 50-1012 | 12 Fr x 21" | Υ | 7 | Angled tip; Open; With control valve; Without tubing; Straight, Nonlatex; Sterile |
| SOT-0055 | Suction Catheter 14Fr | Med-Rx | CHS 50-1014 | 14 Fr x 21" | Υ | 7 | |
| SOT-0057 | Yankauer Suction Bulb Tip | Medline | DYND50130 | Each | Υ | 1 | |

TRACHEOSTOMY PRODUCTS

NOTE: Additional products are available in the TRACHEOTOMY / DRAIN SPONGES category in the Dressings, Bandages, Tape, Compression & Wraps section; Non fenestrated inner cannulas (4IC65; 5IC70; 6IC75; 7IC80; 8IC85) are NOT compatible with fenestrated trach tubes.

| SOT-0032 | Tracheostomy Care Kit | Med-Rx | 85-4037 | Each | Υ | | Includes Disposable Tray, 4x4 Gauze, Plastic transfer Forcep, Pipe Cleaners, Cotton Tipped Applicators, Trach Sponges, Trach Brush, Twill Tape |
|----------|--|-----------------|--------------------|-------------------------------|---|----|--|
| SOT-0034 | Tracheostomy Ties One Size | Shiley | Medtronic TTH | Each | Υ | 1 | Fits both adults and pediatric patients with neck size up to 18" |
| SOT-0014 | Trach Inner Cannula Disposable Beveled Tip for 6.5mm Tube Green | Shiley | Medtronic 4IC65 | I.D. x Length 5.5mm x 62mm | Υ | 71 | For use with Tracheostomy Tubes Green 4CN65H or 4UN65H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use |
| SOT-0102 | Trach Inner Cannula Disposable Beveled Tip for 7.0mm Tube Orange | Shiley Flexible | Medtronic 5IC70 | I.D. x Length 6.0mm x 68mm | Υ | 71 | For use with Tracheostomy Tube Orange 5CN70H or 5UN70H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use |
| SOT-0015 | Trach Inner Cannula Disposable Beveled Tip for 7.5mm Tube Light Blue | Shiley Flexible | Medtronic 6IC75 | I.D. x Length 6.5mm x 74mm | Υ | | For use with Tracheostomy Tube Light Blue 6CN75H or 6UN75H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use |
| SOT-0103 | Trach Inner Cannula Disposable Beveled Tip for 8.0mm Tube Purple | Shiley Flexible | Medtronic 7IC80 | I.D. x Length 7.0mm x 77mm | Υ | | For use with Tracheostomy Tube Purple 7CN80H or 7UN80H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use |
| SOT-0018 | Trach Inner Cannula Disposable Beveled Tip for 8.5mm Tube Navy Blue | Shiley Flexible | Medtronic 8IC85 | I.D. x Length 7.5mm x 79mm | Υ | | For use with Tracheostomy Tube Navy 8CN85H or 8UN85H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | АРО | Weekly Maximum | Notes | |
|---------------------|--|------------------|--|-----------------|--------|-------------------|--|--|
| SOT-0080 | Fenestrated Inner Cannula Disposable Size 4 | Shiley Legacy | Medtronic 4DICFEN | 4.9mm x 62mm | Υ | 21 | For use with Tracheostomy Tubes 4DFEN, 4DCFN | |
| SOT-0081 | Fenestrated Inner Cannula Disposable Size 6 | Shiley Legacy | Medtronic 6DICFEN | 6.0mm x 62mm | Υ | 21 | For use with Tracheostomy Tubes 6DFEN, 6DCFN | |
| SOT-0082 | Fenestrated Inner Cannula Disposable Size 8 | Shiley Legacy | Medtronic 8DICFEN | 7.5mm x 79mm | Υ | 21 | For use with Tracheostomy Tubes 8DFEN, 8DCFN | |
| The following | g Tracheostomy tubes are to be ordered only w | hen suctioning | is required for FE | ENESTRATED trac | heosto | my tubes in | situ | |
| 15(1)-(1)(1// | Inner Cannula Disposable Size 4 for fenestrated tracheostomy tubes | Shiley Legacy | Medtronic 4DIC | 4.9mm x 62mm | Υ | <i>-</i> | For use with fenestrated Tracheostomy Tubes 4DICFEN for suctioning only | |
| 15(1)1-(1)1/X | Inner Cannula Disposable Size 6 for fenestrated tracheostomy tubes | Shiley Legacy | Medtronic 6DIC | 6.0mm x 62mm | Υ | 21 | For use with fenestrated Tracheostomy tubes 6DICFEN for suctioning only | |
| 1 N 1 - 1 1 1 7 M | Inner Cannula Disposable Size 8 for fenestrated tracheostomy tubes | Shiley Legacy | Medtronic 8DIC | 7.5mm x 79mm | Υ | 21 | For use with fenestrated Tracheostomy tubes 8DICFEN for suctioning only | |
| OXYGEN N | MONITORING | | | | | | | |
| EPY-1008 | Oximetry Monitor | | NOTE: Equipment is available for rental as per approved policies from Medigas (Praixair). Please refer to the Medigas Medical Equipment Catalogs for more detailed information before ordering. | | | | | |
| SOT-0099 | Infant Toe Sensor for EPY-1008 O2 Monitor | Masimo Canada | Masimo 2328 | Each | Υ | 1 | Single Patient Use. Includes one adhesive tape. | |
| SOT-0100 | Replacement Tapes for Infant Toe Sensor | Masimo Canada | Masimo 2307 | Each | Υ | 5 | Single Patient Use, Replacement adhesive tape for SOT-0099 | |

PROTECTIVE WEAR

| CHRIS | Description | Brand Name | Manufacturer | Sizo | ΛDO | Weekly | Notes |
|------------|-------------|---------------|--------------|------|-----|---------|-------|
| Order Code | Description | Dialiu Naille | Code | Size | AFO | Maximum | Notes |

PROTECTIVE WEAR: Personal protective equipment (PPE) is authorized for Home and Community Care Support Services North East staff to protect employees from exposure to workplace hazards as per policy. Service Providers are expected to provide their own personal protective equipment to their own staff to protect from infectious diseases. Service Providers may be provided PPE if handling hazarous chemicals/drugs for approved patients. Non-Sterile PPE may be authorized for contracted Service Providers only when the on-going patient-specific treatment requires protective wear be worn to minimize risk to patient as approved by the Care Coordinator. Clean PPE is standard in community care, and sterile PPE supplies are only to be authorized when treatments require a sterile technique as per approved policies or clinical guidelines.

| SOT-0004 | Disposable Isolation Gown Regular/Large | Medline | Medline CRI4000 | Each | N | 1 per Visit | Single use disposable; For use for protection from chemotherapy or other cytotoxic medications in liquid | |
|----------|---|-------------------|----------------------------|--------------------|---|-------------|--|--|
| SOT-0098 | Disposable Chemo Gown Regular/Large | Medline | Medline NON27457 | Each | N | 1 per Visit | form (Intravenous, subcutaneous, intramuscular, liquid oral, etc.). NOTE: This face mask is not appropriate if | |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | Primaguard 120 | priMED Medical PG4-1092 | Each | N | 1 per Visit | the cytotoxic medication is aerosolized - Use a fit-tested N95 mask provided by the nurse's employer. | |
| SOT-0019 | Eye Protection Safety Glasses with Side Shields | Medical Mart | 320-NON24777V | Each | N | 1 | For procedures and patient-care activities likely to generate splashes or sprays of blood or body fluids or cytotoxic medications in liquid form | |
| SOT-0020 | Sterile Nitrile Gloves Small (Pair) | Flexam | Cardinal Health N8830 | Each Pair | Υ | 3 | Non-latex, powder free, nitrile examination gloves with 11 in. cuff. Approved for use with cytotoxic | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | Flexam | Cardinal Health N8831 | Each Pair | Υ | 3 | medications. Nursing providers are requried to provide their staff with | |
| SOT-0022 | Sterile Nitrile Gloves Large (Pair) | Flexam | Cardinal Health N8832 | Each Pair | Υ | 3 | non-sterile gloves. Sterile gloves will only to be provided to nursing service provider personnel when a sterile | |
| SOT-0023 | Sterile Nitrile Gloves Extra Large (Pair) | Flexam | Cardinal Health N8833 | Each Pair | Υ | 3 | technique is required for patient care as approved by the Care Coordinator | |
| SOT-0012 | Sterile Nitrile Gloves Medium (Box) | Flexam | Cardinal Health N8831 | Box of 40 Pairs | N | 1 | Non-latex, powder free, nitrile examination gloves Bulk orders available Clinics and Residental Hospice only. Approved for use with cytotoxic medications. | |
| SOT-0093 | Non-Sterile Nitrile Gloves Small (Pair) | Esteem XP | Cardinal Health N8851XP | Each Pair | N | 7 | Non-latex, powder free, nitrile examination gloves. ASTM F 739 Chemicals tested, but not for use with | |
| SOT-0094 | Non-Sterile Nitrile Gloves Medium (Pair) | Esteem XP | Cardinal Health N8852XP | Each Pair | N | 7 | chemotherapy. Nursing providers are required to provide their staff with | |
| SOT-0096 | Non-Sterile Nitrile Gloves Large (Pair) | Esteem XP | Cardinal Health N8853XP | Each Pair | N | 7 | non-sterile gloves. | |
| SOT-0095 | Non-Sterile Nitrile Gloves Medium (Box) | Esteem XP | Cardinal Health N8852XP | Box of 50 | N | 1 | Bulk orders only available for Home and Community Care North East clinical staff or Clinics or Hospices | |

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests ALGOMA SUB-REGION (Questions or issues call: 705-949-1808 / 1-800-668-7705 Ext: 3269)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|--|---|--|--|
| URBAN | Soo Medical | 656 Great Northern Road Sault Ste. Marie, ON P6B 4Z9 705-574-0228 | Medicine Shoppe's Infusion (IV) Equipment & Supplies | |
| | | Fax: 705-574-0231 | Robinson's Medical Supplies | |
| URBAN | The Medicine Shoppe | 316 Second Line West Sault Ste. Marie, ON P6C 2J5 705-253-0720/877-374-6938 | Medicine Shoppe's Infusion (IV) Equipment & Supplies | |
| | D 14 100 1 | Fax: 705-256-7656 | Robinson's Medical Supplies | |
| URBAN | ParaMed Clinic | 477 Queen Street East Suite 100 Sault Ste. Marie, ON P6A 1Z5 705-759-0110 | Medicine Shoppe's Infusion (IV) Equipment & Supplies | For ParaMed Clinic Patients only |
| | | | Robinson's Medical Supplies | *Not a public depot |
| URBAN | Algoma Residential Hospice Care (ARCH) | 229 Fourth Line West Sault Ste. Marie, ON P6A 0B5 705-942-1556 | Medicine Shoppe's Infusion (IV) Equipment & Supplies | For Hospice Residents only *Not a public depot |
| DUDAL | | 4400 5: 1 1 0: 1 | Robinson's Medical Supplies | |
| RURAL | Matthews Memorial Hospital | 1180 Richards Street Richards Landing, ON P0R 1J0 705-246-2570 | Robinson's Medical Supplies | |
| RURAL | Main Street Pharmacy | 156 Main Street Thessalon, ON P0R 1L0 705 842 2322 | Medicine Shoppe's Infusion (IV) Equipment & Supplies | |
| | | | Robinson's Medical Supplies | |
| RURAL | Pit Stop Gas Bar & Confectionary | 3245 Highway 17 East Echo Bay, ON P0S 1C0 705-248-3333 | Robinson's Medical Supplies | |
| RURAL | Goulais River Country Store | 2771 Highway 17 North Goulais River, ON P0S 1E0 705-649-3811 | Robinson's Medical Supplies | |
| RURAL | Elliot Lake PharmaChoice | 267 Hillside Drive South Elliot Lake, ON P5A 1N7 705-848-0004 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Home Hardware | 275 Imperial Street Massey, ON P0P 1P0 705-865-2135 | Robinson's Medical Supplies | |
| RURAL | Fenlon's Pharmacy | 88 Broadway Avenue Wawa, ON POS 1K0 705-856-2800 | Robinson's Medical Supplies | |
| RURAL | Patient Service Treatment Location - Hornepayne | Patient Service Treatment Location | Medicine Shoppe's Infusion (IV) Equipment & Supplies Robinson's Medical Supplies | Patients living in the pre- defined Hornepayne area are pre-authorized for home delivery as there is no available depot. |
| RURAL | Patient Service Treatment Location - Blind River | Patient Service Treatment Location | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | Patients living in the predefined Blind River area are pre-authorized for home delivery because there is no available depot. |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests NIPISSING-TEMISKAMING SUB-REGION (Questions or issues call: Robinson's Pharmasave 705-495-4515 - Option 3. NOTE: For Sturgeon Falls medical supplies only call 705-566-5551 - Option 4)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|--|--|--|--|
| URBAN | VON Ambulatory Nursing Clinic | 150 First Avenue West North Bay, ON P1N 3B9 | Robinson's Medical Supplies & Infusion (IV) Equipment & | For VON Clinic Patients only |
| | | 705-472-8050 Fax: 705-495-6691 | Supplies | *Not a public depot |
| URBAN | Robinson's Pharmasave | 2547 Trout Lake Road North Bay, ON P1B 7S8 705-495-4515 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| URBAN | CarePartners Ambulatory Nursing Clinic | 195 First Avenue West North Bay, ON P1B 3B8 1-866-886-4172 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | For CarePartners Clinic Patients only |
| URBAN | Nipissing Serenity Hospice | 799 John Street North Bay, ON P1B 8T2 | Robinson's Medical Supplies & Infusion (IV) Equipment & | *Not a public depot For Hospice Residents Only |
| RURAL | Pharmacie Nipissing Pharmacy, PharmaChoice | 705-995-3377 26 Principale Street East Verner, ON P0H 2M0 705-581-2013 | Supplies Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | *Not a public depot |
| RURAL | Mattawa IDA | 300 Main Street Mattawa, ON P0H 1V0 705-744-5981 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | South River Pharmacy | 311 Hwy 124 South River, ON P8A 1X0 705-386-0636 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Loblaws Pharmacy (No Frills) | 12035 Hwy 17 Sturgeon Falls, ON P2B 2S6 705-753-5850 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Clement's IDA | 497 Main Street Powassan, ON P0H 1Z0 705-724-3523 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Argyle Nursing Station | 11851 Hwy 522 Port Loring, ON P0H 1Y0 705-757-1717 | Robinson's Medical Supplies | |
| RURAL | Arnica Medical Pharmacy | 162 Huston Street Burks Falls, ON P0A 1C0 705-382-2111 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Sundridge Guardian Pharmacy | 107 Main Street Sunridge, ON P0A 1Z0 705-384-5335 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Closing the Gap Healthcare Clinic | 395 Centre Street North Suite 202 Huntsville, ON P1H 1H2 | Robinson's Medical Supplies | For Closing the Gap Clinic patients only |
| | | 705-788-9355 | | *Not a public depot |

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

NIPISSING-TEMISKAMING SUB-REGION (Questions or issues call Kirkland Lake office at 1- 888-533-2222 Ext: 3269 or call

Robinson's Pharmasave 705-495-4515 (Option 3)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|--------------------------------|--|--|---|
| RURAL | Pharmasave BDR Drug Mart | 15 Station Road North Kirkland Lake, ON P2N 3J1 705-567-3315 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | J&S Variety Store | 91 Pine Street Elk Lake, ON P0J 1G0 705-678-2246 | Robinson's Medical Supplies | |
| RURAL | Elk Lake Nursing Station | 71 Front Street Elk Lake, ON P0J 1G0 705-678-2215 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Larder Lake Co-Op | 25 Godfrey Street Larder Lake, ON P0K 1L0 705-643-2327 | Robinson's Medical Supplies | |
| RURAL | Larder Lake Pharmasave | 15 Godfrey Street Larder Lake, ON P0K 1L0 705-643-2400 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Dogstar Service | 418 Moyneur Avenue Matachewan, ON P0K 1M0 705-565-2262 | Robinson's Medical Supplies | *IV Meds and supplies delivered to patient's treatment location |
| RURAL | CarePartners Nursing Clinic | 883361 Highway 65 North New Liskeard, ON P0J 1P0 1-866-886-4172 ext:5351 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | For CarePartners' Clinic Patients only *Not a public depot |
| RURAL | Loblaw Pharmacy | 55 Scott Street New Liskeard, ON P0J 1P0 705-647-8728 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | Trock a passio dopor |
| RURAL | Latchford Town Office | 10 Main Street Latchford, ON P0J 1N0 705-676-2416 | Robinson's Medical Supplies | |
| RURAL | Temagami Family Health Team | 17 O'Connor Drive Temagami, ON P0H 2H0 705-569-3244 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Marshall's Pharamsave | 43 Third Street Englehart, ON P0J 1H0 705-544-7878 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Cote's Variety Store | 16 Tenth Street East Earlton, ON P0J 1E0 705-563-2022 | Robinson's Medical Supplies | |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

COCHRANE SUB-REGION (Questions or issues call: 1-705-269-0001)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|--|--|--|--|
| URBAN | Bayshore Home Health Clinic | 202-119 Pine Street Timmins, ON P4N 2K3 705-268-6088 / 866-227-3661 Fax: 705-268-9085 | Bayshore Rx IV Supplies & Equipment Robinson's Medical Supplies | Bayshore Clinic patient's |
| RURAL | Iroquois Falls Pharmacy | Iroquois Falls Pharmacy Rx Drug Mart 201 Ambridge Drive Iroquois Falls, ON P0G 1G0 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment are delivered to patient's treatment location * Public depot for Robinson's Medical Supplies only |
| RURAL | Matheson Willis Pharmacy | Willis Pharmacy 414 5th Avenue Matheson, ON P0K 1N0 705-273-2202 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment are delivered to patient's treatment location * Public depot for Robinson's Medical Supplies only |
| RURAL | Smooth Rock Falls Pharmacy | 133 5th Street Smooth Rock Falls, ON P0L 2B0 705-338-4470 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment is delivered to patient's treatment location * Public depot for Robinson's Medical Supplies only |
| RURAL | Cochrane Pharmasave | Cochrane Pharmasave 134A Third Street Cochrane, ON P0L 1C0 705-272-1333 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment are delivered to patient's treatment location * Public depot for Robinson's Medical Supplies only |
| RURAL | Hearst | Pharmacie Novena 1509 Ontario 11 West Hearst, ON P0L 1N0 705-372-1212 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment are delivered to patient's treatment location * Public depot for Robinson's Medical Supplies only |
| RURAL | Kapuskasing | Kapuskasing Family Health Team 142C Progress Crescent Kapuskasing, ON P5N 3H6 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment are delivered to patient's treatment location * Public depot for Robinson's Medical Supplies only |
| RURAL | Missinaibi-Mattice Clinic | 500 Hwy 11 East Mattice, ON P0L 1T0 705-364-2161 | Robinson's Medical Supplies | All Bayshore Rx IV Supplies & Equipment is delivered to patient's home as there is no |
| RURAL | Ontario Health atHome Moosonee Office | PO Box 730, 48 Revillion Street Moosonee ON, P0L 1Y0 705-336-2233 Fax: 705-336-2919 | Robinson's Medical Supplies | available public depot. * Public depot for Robinson's Medical Supplies only |
| RURAL | Cavalier Motel | 98 Monk Street Chapleau, ON P0M 1K0 705-864-0841 | Robinson's Medical Supplies | |
| RURAL | Jr Corner Store Gogama | 41 Poupore Street Gogama, ON, P0M 1W0 705-894-2540 Attention: Gogama Nursing Clinic | Robinson's Medical Supplies | |

Private-Pay Patient Deliveries: Robinson's to coordinate for Patients when requested

Note: All IV orders from Bayshore Rx are delivered directly to the patient's treatment location

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

SUDBURY-MANITOULIN-PARRY SOUND SUB-REGION

(Questions or issues call: Robinson's Pharmasave (705-566-5551 Option 4/1-800-598-2988 Option 4)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|---------------------------------------|--|--|---|
| URBAN | Health Care Pharmacy Pharmasave | 1276 Lasalle Boulevard Sudbury, ON P3A 1Y8 705-566-5592 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| URBAN | Bayshore Ambulatory Nursing Clinic | 2120 Regent Street, Unit 8 Sudbury, ON P3E 5S8 705-523-6668 | Robinson's Medical Supplies and IV Equipment & IV Supplies | For Bayshore Clinic Patients only |
| URBAN | Larch St. Place | 200 Larch Street Sudbury, ON P3E 1C5 705-673-3721 | Robinson's Medical Supplies | For Sudbury District Nurse Practitioner Clinic patients only |
| URBAN | Daffodil Lodge | 41 Ramsey Lake Road Sudbury, ON P3E 5J1 705-523-7333 | Robinson's Medical Supplies | For patients residing at Daffodil Lodge only |
| URBAN | ICAN (Haig) | 582 Haig Street Sudbury, ON P3C 1C4 | Robinson's Medical Supplies | For ICAN patients only |
| URBAN | ICAN (Mont Adam) | 117 Mont Adam Sudbury, ON P3B 1E1 | Robinson's Medical Supplies | For ICAN patients only |
| URBAN | Sudbury Maison McCulloch Hospice | 1028 South Bay Sudbury, ON P3E 6J7 705-674-9252 | Robinson's Medical Supplies and IV Equipment & IV Supplies | For Hospice Residents only |
| URBAN | Health Sciences North | 41 Ramsey Lake Road Sudbury, ON P3E 5J1 | Robinson's Medical Supplies and IV Equipment & IV Supplies | For Ontario Health atHome Stock only |
| | | | | Not a depot for Patient pick-up |
| URBAN | Chelmsford Guardian Pharmacy | 46 Main Street West Chelmsford, ON P0M 1L0 705-855-2371 | Robinson's Medical Supplies | |
| URBAN | Uptown Cleaners | 1610 Regent Street Sudbury, ON P3E 705-522-0400 | Robinson's Medical Supplies | |
| URBAN | Coniston Guardian | 24 Amanda Street Coniston, ON P0M 1M0 705-694-5346 | Robinson's Medical Supplies | |
| URBAN | Wilson Pharmacy Pharmasave | 10 Serpentine Street Copper Cliff, ON P0M 1N0 705- 682-2841 | Robinson's Medical Supplies | |
| URBAN | Nickel Centre Pharmacy, Pharmasave | 3098 Falconbridge Highway Garson, ON P3L 1P5 705-693-2738 | Robinson's Medical Supplies | |
| URBAN | Val Est Pharmacy, Pharmasave | 3179 Highway 69 North Unit # 1 Val Caron, ON P3N 1S5 705-897-1867 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| URBAN | Walden Medicine Shoppe | 135 Regional Road 24 Walden Plaza Lively, ON P3Y 1J1 705-692-5868 | Robinson's Medical Supplies | |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

SUDBURY-MANITOULIN-PARRY SOUND SUB-REGION

(Questions or issues call: Robinson's Pharmasave (705-566-5551 Option 4/1-800-598-2988 Option 4)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|--|---|--|-------|
| URBAN | Capreol Rx | 21 Young Street Capreol, ON P0M 1H0 705-858-8888 Fax: 705-858-9888 | Robinson's Medical Supplies | |
| URBAN | Dowling Pharmacy | 30 Main Street Dowling, ON 705-855-4562 | Robinson's Medical Supplies | |
| URBAN | Levack Pharmacy | 2 – 1 st Avenue North Levack, ON P0M 2C0 705-966-3612 | Robinson's Medical Supplies | |
| RURAL | Hagar General Store | 7147 Highway 17 East Hagar, ON 705-599-2010 | Robinson's Medical Supplies | |
| RURAL | Noelville Pharmacy | 27 Street David North Noelville, ON P0O 2N0 705-898-2999 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | St. Charles Cooperative Store | 25 Main Street Charles, ON P0M 2W0 705-867-2000 | Robinson's Medical Supplies | |
| RURAL | Pitfield General Store | 7 Channel Killarney, ON P0M 2A0 705-287-2872 | Robinson's Medical Supplies | |
| RURAL | Robinsons Pharmasave | 119 Tudhope Espanola, ON P5E 1S6 705-869-1561 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Home Hardware | 275 Imperial Street Massey, ON P0P 1P0 705-865-2135 | Robinson's Medical Supplies | |
| RURAL | Little Current Guardian | 41 Meredith Street Little Current, ON P0P 1K0 705-368- 2530 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Little Current Hospital Nursing Station | 11 Meredith Street West Little Current, ON P0P 1K0 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Manitowaning Pharmacy | 147 Arthur Street Manitowaning, ON P0P 1N0 705-859-2323 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Mindemoya Pharmacy | Guardian Pharmacy 11 Yonge Street Mindemoya, ON P0P 1S0 705-377-5411 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Robertson IDA | 39 Meredith Street Gore Bay, ON P0P 1H0 705-282-2147 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Central Guardian Pharmacy | 45 Meredith Street Gore Bay, ON P0P 1H0 705-282-2409 | Robinson's Medical Supplies | |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

SUDBURY-MANITOULIN-PARRY SOUND SUB-REGION

(Questions or issues call: Robinson's Pharmasave (705-566-5551 Option 4/1-800-598-2988 Option 4)

| Service Area | Depot Name | Location | Service Type | Notes |
|-----------------|--|--|--|--|
| RURAL | Britt Nursing Station | Box 75 991 Riverside Drive Britt, ON P0G 1A0 705-383-2375 | Robinson's Medical Supplies | |
| RURAL | Duck Rock General Store | 2093 124 Highway Dunchurch, ON P0A 1G0 705-389-3322 | Robinson's Medical Supplies | |
| RURAL | Lane Family Pharmasave | 4 Pine Drive Parry Sound, ON P2A 3B8 705-746-2187 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Rosseau General Store | 1 Rice St. Rosseau, ON P0C 1J0 705-732-4479 | Robinson's Medical Supplies | |
| RURAL | Ontario Health atHome West Parry Sound Hospital | 6 Albert Street Parry Sound, ON P2A 3A4 | Robinson's Medical Supplies | For Ontario Health atHome Stock only Not a depot for Patient pick-up |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life



SUDBURY OFFICE
40 Elm Street, Suite 41-C
Sudbury ON P3C 1S8
Tel : 705 522-3461
Toll Free : 1 800 461-2919
Fax : 705 522-3855
TTY: 711 (ask operator 1-888-533-2222)

MARKA pelhin on ca

BUREAU DE SUDBURY
40, rue Elm, bureau 41-C
Sudbury ON P3C 1S8
Téléphone: 705 522-3461
Téléphone: 705 522-3461
Telécopieur : 705 522-3855
TTY: 711 (ask operator 1-888-533-2222)

MARKA pelhin on ca

WWW pelhin on ca

www.nelhin.on.ca

www.nelhin.on.ca